

Product datasheet for **RC216002L4V**

SLC25A33 (NM_032315) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC25A33 (NM_032315) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC25A33
Synonyms:	BMSC-MCP; PNC1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032315
ORF Size:	963 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216002).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032315.1
RefSeq Size:	1476 bp
RefSeq ORF:	966 bp
Locus ID:	84275
UniProt ID:	Q9BSK2
Cytogenetics:	1p36.22
Domains:	mito_carr
Protein Families:	Transmembrane



[View online »](#)

MW: 35.4 kDa

Gene Summary: SLC25A33 belongs to the SLC25 family of mitochondrial carrier proteins (Haitina et al., 2006 [PubMed 16949250]).[supplied by OMIM, Mar 2008]