

## Product datasheet for **RC216039L4V**

### IRF5 (NM\_001098631) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	IRF5 (NM_001098631) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IRF5
Synonyms:	interferon regulatory factor 5; SLEB10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001098631
ORF Size:	1464 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216039).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001098631.1</a> , <a href="#">NP_001092101.1</a>
RefSeq Size:	2741 bp
RefSeq ORF:	1466 bp
Locus ID:	3663
Cytogenetics:	7q32.1
Protein Families:	Transcription Factors
Protein Pathways:	Toll-like receptor signaling pathway
MW:	54.8 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the interferon regulatory factor (IRF) family, a group of transcription factors with diverse roles, including virus-mediated activation of interferon, and modulation of cell growth, differentiation, apoptosis, and immune system activity. Members of the IRF family are characterized by a conserved N-terminal DNA-binding domain containing tryptophan (W) repeats. Alternative promoter use and alternative splicing result in multiple transcript variants, and a 30-nt indel polymorphism (SNP rs60344245) can result in loss of a 10-aa segment. [provided by RefSeq, Dec 2016]