

## Product datasheet for **RC216058**

### TCF7L2 (NM\_030756) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	TCF7L2 (NM_030756) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TCF7L2
Synonyms:	TCF-4; TCF4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RC216058 representing NM\_030756  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGCCGACGTGAACGGCGGTGGAGGGGATGACCTAGGCGCCAACGACGAACTGATTTCTTCAAAGACG  
 AGGGCGAACAGGAGGAGAAGAGCTCCGAAAACCTCTCGGCAGAGAGGGATTTAGCTGATGTCAAATCGTC  
 TCTAGTCAATGAATCAGAAACGAATCAAAAACAGCTCCTCCGATTCGAGGCGGAAAAGACGGCCTCCGCC  
 CGCTCCGAAAAGTTCCGAGACAAATCCCGGAAAGTTTGAAGAAGCGGCCAAGAGGCAAGATGGAGGGC  
 TCTTTAAGGGGCCACCGTATCCCGGCTACCCCTTCATCATGATCCCCGACCTGACGAGCCCTACCTCCC  
 CAACGGATCGCTCTCGCCACCGCCCGAACCTATCTCCAGATGAAATGGCCACTGCTTGTGTCCAGGCA  
 GGGAGCCTCCAGAGTAGACAAGCCCTCAAGGATGCCCGGTCCCCATACCGGCACACATTGTCTTAACA  
 AAGTGCCAGTGGTGCAGCACCCCTCACCATGTCCACCCCTCACGCCTTTATCACGTACAGCAATGAACA  
 CTTACGCGGGAAACCCACCTCCACACTTACCAGCCGACGTAGACCCAAAACAGGAATCCACGCGCT  
 CCGCACCTCCAGATATATCCCGTATTACCCACTATCGCTGGCACCGTAGGACAAATCCCCATCCGC  
 TAGGATGGTTAGTACCACAGCAAGGTCAACCAGTGTACCCAATCACGACAGGAGGATTCAGACACCCCTA  
 CCCCACAGCTCTGACCGTCAATGCTTCCATGTCCAGGTTCCCTCCCCATATGGTCCCACCACATCATACG  
 CTACACACGACGGGCAATCCGCATCCGGCCATAGTCACACCAACAGTCAAACAGGAATCGTCCCAGAGTG  
 ATGTCGGCTCACTCCATAGTTCAAAGCATCAGGACTCCAAAAGGAAGAAGAAAAGAAAGAGCCCAT  
 AAAGAAACCTCTTAATGCATTATGTTGTATATGAAGGAAATGAGAGCAAAGGTCGTAGCTGAGTGCAGC  
 TTGAAAGAAAGCGCGCCATCAACCAGATCCTTGGCGGAGGTGGCATGCACTGTCCAGAGAAGAGCAAG  
 CGAAATACTACGAGCTGGCCCGGAAGGAGCGACAGCTTCATATGCAACTGTACCCCGGCTGGTCCGCGG  
 GGATAACTATGGAAAGAAGAAGAAGAGGAAAGGACAAGCAGCGGAGAGACCAATGAACACAGCGAA  
 TGTTTCCTAAATCCTTGCCTTCTACTTCTCCGATTACAGACCTGAGCGCTCCTAAGAAATGCCAGCGC  
 GCTTTGGCCTTGATCAACAGAATAACTGGTGGCGCCCTTGAGGAGAAAAAAAAGTGCCTTGCCTACAT  
 ACAAGGTGAAGGCAGCTGCCTCAGCCCAACCTCTTCAGATGGAAGCTTACTAGATTGCCTCCCCCTCC  
 CCGAACCTGCTAGGCTCCCTCCCGAGACGCCAAGTACAGACTGAGCAGACCCAGCCTCTGTCTGTGT  
 CCCTGAAGCCGACCCCTGGCCACCTGTCCATGATGCCTCCGCCACCCGCCCTCTGTCTGTGAGGC  
 CACCCACAAGGCCTCCGCCCTGTGCCAACGGGGCCCTGGACCTGCCCCAGCCGCTTTCAGCCTGCC  
 GCCCCTCCTCATCAATTGCACAGCCGTGACTTCTTCTTACATTCCCACAGCTCCTGGCCGGGACCC  
 AGCCCCAGCCGCTGTCGCTCGTACCAAGTCTTTAGAA

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**

>RC216058 representing NM\_030756  
 Red=Cloning site Green=Tags(s)

MPQLNGGGDDLGANDELSFKDEGEQEEKSSSENSAERDLADVKSSLVNESETNQSSSDSEARRPPP  
 RSESFDRKSRESLEEAQRQDGLFKGPPYPGYPFIMIPDLTSPYLPNGSLSPARTYTLQMKWPLLDVQA  
 GSLQSRQALKDARSPSPAHIVSNKVPVQHPHHVHPLTPLITYSNEHFTPGNPPPHLPADVDPKTIIPRP  
 PHPPDISPYYPSPGTVGQIPHPLGWLVPQQGQPVYPIITGGFRHPYPTALTVNASMSRFPHPMVPPHHT  
 LHHTGIPHPAIVTPTVKQESSQSDVGLSHSSKHQDSKKEEKKKPHIKKPLNAFMLYMKEMRAKVVAECT  
 LKESAAINQILGRRWHALSREEQAKYYELARKERQLHMQLYPGWSARDNYGKKKKRKRDKQPGETNEHSE  
 CFLNPCLSLPPIIDLAPKKCRARFGLDQNNWCGPCRKKKCVRYIQEGSCLSPSSDGLLDSPPPS  
 PNLLGSPPRDAKSQTEQTQPLSLKPDPLAHL SMMPPPALLLAEATHKASALCPNGALDLPPAALQPA  
 APSSSIAQPSTSSLHSHSSLAGTQPQPLSLVTKSLE


SGP**TRRRLEQKLI**SEEDLAANDILDYKDDDDKV

**Chromatograms:**

[https://cdn.origene.com/chromatograms/mg4242\\_g08.zip](https://cdn.origene.com/chromatograms/mg4242_g08.zip)

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM\_030756

ORF Size: 1788 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_030756.5](#)

RefSeq Size: 2444 bp

RefSeq ORF: 1791 bp

Locus ID: 6934

**Cytogenetics:** 10q25.2-q25.3

**Domains:** HMG


**Protein Families:** Druggable Genome, Transcription Factors

**Protein Pathways:** Acute myeloid leukemia, Adherens junction, Arrhythmogenic right ventricular cardiomyopathy (ARVC), Basal cell carcinoma, Colorectal cancer, Endometrial cancer, Melanogenesis, Pathways in cancer, Prostate cancer, Thyroid cancer, Wnt signaling pathway


**MW:** 65.1 kDa

**Gene Summary:** This gene encodes a high mobility group (HMG) box-containing transcription factor that plays a key role in the Wnt signaling pathway. The protein has been implicated in blood glucose homeostasis. Genetic variants of this gene are associated with increased risk of type 2 diabetes. Several transcript variants encoding multiple different isoforms have been found for this gene.[provided by RefSeq, Oct 2010]

**Product images:**



Circular map for RC216058



Western blot validation of overexpression lysate (Cat# [LY403074]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC216058 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).