

Product datasheet for RC216058L2

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

TCF7L2 (NM_030756) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: TCF7L2 (NM_030756) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: TCF7L2

Synonyms: TCF-4; TCF4

Mammalian Cell None

Selection:

Vector: pLenti-C-mGFP (PS100071)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide


The ORF insert of this clone is exactly the same as(RC216058).

Sequence:

Restriction Sites: Sgfl-Rsrll

Cloning Scheme:





 $[\]ensuremath{^*}$ The last codon before the Stop codon of the ORF.

ACCN: NM_030756

ORF Size: 1788 bp





TCF7L2 (NM_030756) Human Tagged Lenti ORF Clone - RC216058L2

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 030756.1</u>

RefSeq Size: 2444 bp RefSeq ORF: 1791 bp Locus ID: 6934

Cytogenetics: 10q25.2-q25.3

Domains: HMG

Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Acute myeloid leukemia, Adherens junction, Arrhythmogenic right ventricular

cardiomyopathy (ARVC), Basal cell carcinoma, Colorectal cancer, Endometrial cancer,

Melanogenesis, Pathways in cancer, Prostate cancer, Thyroid cancer, Wnt signaling pathway

MW: 65.1 kDa

Gene Summary: This gene encodes a high mobility group (HMG) box-containing transcription factor that plays


a key role in the Wnt signaling pathway. The protein has been implicated in blood glucose homeostasis. Genetic variants of this gene are associated with increased risk of type 2

diabetes. Several transcript variants encoding multiple different isoforms have been found for

this gene.[provided by RefSeq, Oct 2010]



Product images:



Circular map for RC216058L2