

Product datasheet for **RC216082A1V**

Human CTNNBIP1 (NM_020248) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human CTNNBIP1 (NM_020248) AAV Particle

Tag: Myc-DDK

Symbol: CTNNBIP1

Synonyms: ICAT

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC216082 representing NM_020248
Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAACCGCGAGGGAGCTCCCGGAAGAGTCCGGAGGAGATGTACATTCAGCAGAAGGTCCGAGTGCTGC
TCATGCTGCGGAAGATGGGATCAAACCTGACAGCCAGCGAGGAGGAGTTCCTGCGCACCTATGCAGGGGT
GGTCAACAGCCAGCTCAGCCAGCTGCCTCCGCACTCCATCGACCAGGGTGCAGAGGACGTGGTATGGCC
TTTTCCAGGTCGGAGACGGAAGACCGGAGGCAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC216082 representing NM_020248
Red=Cloning site Green=Tags(s)

MNREGAPGKSPEEMYIQQKVRVLLMLRKMGSNLTASEEEFLRTYAGVVNSQLSQLPPHSIDQGAEDVVMA
FSRSETEDRRQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2

ACCN: NM_020248

ORF Size: 243 bp

Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_020248.2</u>
RefSeq Size:	2995 bp
RefSeq ORF:	246 bp
Locus ID:	56998
UniProt ID:	<u>Q9NSA3</u>
Cytogenetics:	1p36.22
MW:	9 kDa