

Product datasheet for RC216106L2

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

GLB1 (NM_000404) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: GLB1 (NM_000404) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: GLB1

Synonyms: EBP; ELNR1; MPS4B

Mammalian Cell None

Selection:

Vector: pLenti-C-mGFP (PS100071)

E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC216106).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_000404

ORF Size: 2030 bp

GLB1 (NM_000404) Human Tagged Lenti ORF Clone - RC216106L2

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 000404.1</u>

RefSeq Size:2409 bpRefSeq ORF:2034 bpLocus ID:2720

Cytogenetics: 3p22.3

UniProt ID:

Domains: Glyco_hydro_35

Protein Families: Druggable Genome

Protein Pathways: Galactose metabolism, Glycosaminoglycan degradation, Glycosphingolipid biosynthesis -

ganglio series, Lysosome, Metabolic pathways, Other glycan degradation, Sphingolipid

metabolism

P16278

MW: 75.9 kDa


Gene Summary: This gene encodes a member of the glycosyl hydrolase 35 family of proteins. Alternative

splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed to generate the mature lysosomal enzyme. This enzyme catalyzes the hydrolysis of a terminal beta-linked galactose residue from ganglioside substrates and other glycoconjugates. Mutations in this gene may result in GM1-


gangliosidosis and Morquio B syndrome. [provided by RefSeq, Nov 2015]



Product images:



Circular map for RC216106L2



Double digestion of RC216106L2 using Sgfl and Mlul