

Product datasheet for **RC216107L2V**

SLC39A6 (NM_012319) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SLC39A6 (NM_012319) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SLC39A6 |
| Synonyms: | LIV-1; ZIP6 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_012319 |
| ORF Size: | 2265 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC216107). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_012319.3 , NP_036451.3 |
| RefSeq Size: | 3637 bp |
| RefSeq ORF: | 2268 bp |
| Locus ID: | 25800 |
| UniProt ID: | Q13433 |
| Cytogenetics: | 18q12.2 |
| Domains: | Zip |
| Protein Families: | Druggable Genome, Transmembrane |



[View online »](#)

MW: 84.9 kDa

Gene Summary: Zinc is an essential cofactor for hundreds of enzymes. It is involved in protein, nucleic acid, carbohydrate, and lipid metabolism, as well as in the control of gene transcription, growth, development, and differentiation. SLC39A6 belongs to a subfamily of proteins that show structural characteristics of zinc transporters (Taylor and Nicholson, 2003 [PubMed 12659941]).[supplied by OMIM, Mar 2008]