

## Product datasheet for **RC216110L2V**

### GRB 14 (GRB14) (NM\_004490) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	GRB 14 (GRB14) (NM_004490) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GRB 14
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_004490
ORF Size:	1620 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216110).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004490.1</a>
RefSeq Size:	2402 bp
RefSeq ORF:	1623 bp
Locus ID:	2888
UniProt ID:	<a href="#">Q14449</a>
Cytogenetics:	2q24.3
Domains:	RA, SH2, PH
MW:	61 kDa



[View online »](#)

**Gene Summary:**

The product of this gene belongs to a small family of adapter proteins that are known to interact with a number of receptor tyrosine kinases and signaling molecules. This gene encodes a growth factor receptor-binding protein that interacts with insulin receptors and insulin-like growth-factor receptors. This protein likely has an inhibitory effect on receptor tyrosine kinase signaling and, in particular, on insulin receptor signaling. This gene may play a role in signaling pathways that regulate growth and metabolism. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2014]