

## Product datasheet for **RC216111L4V**

### PCGF3 (NM\_006315) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PCGF3 (NM_006315) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PCGF3
Synonyms:	DONG1; RNF3; RNF3A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006315
ORF Size:	726 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216111).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006315.4</a>
RefSeq Size:	5666 bp
RefSeq ORF:	729 bp
Locus ID:	10336
UniProt ID:	<a href="#">Q3KNV8</a>
Cytogenetics:	4p16.3
Domains:	RING
MW:	28.1 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene contains a C3HC4 type RING finger, which is a motif known to be involved in protein-protein interactions. The specific function of this protein has not yet been determined. [provided by RefSeq, Jul 2008]