

Product datasheet for **RC216114**

Mitofilin (IMMT) (NM_001100169) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mitofilin (IMMT) (NM_001100169) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Mitofilin
Synonyms:	HMP; Mic60; MICOS60; MINOS2; P87; P87/89; P89; PIG4; PIG52
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC216114 representing NM_001100169
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTGCGGGCCTGTCTAGTTATCGGGTGTGACCGCCGCCAGAGTTGTCTGTGGGAAGTTTGTCC
 TCCGTCCATTGCGACCATTGCCGACAGATACTCTACTTCAGGCAGCTCTGGGTTGACTACTGGCAAAATTGC
 TGGAGCTGGCCTTTTGTGGTGGAGGATTGGTGGCACTATCCTATATGCCAAATGGGATCCCAT
 TTCCGGGAAAGTGTAGAGAAAACCATACCTTACTCAGACAACTCTTCGAGATGGTTCTGGTCTGCAG
 CTTATAATGTTCCATTGCCAAAGAAATCGATTTCAGTCCGGTCCACTAAAAATCTCTAGTGTATCAGAAGT
 AATGAAAGAACTAAACAGCCTGCCTCACAACCTCAAAAAACAAAAGGGAGATACTCCAGCTTCAGCAACA
 GCACCTACAGAAGCGGCTCAAATATTTCTGCAGCAGGTGATACCCTGTCGGTCCAGCCCTGCAGTTC
 AGCCTGAGGAATCTTAAAAACTGATCACCTGAAATTTGGTGAAGGAAAACCCACACCTGCACCTTCAGA
 AGCATCCTCATCTCTATAAGGGAGCGACCACCTGAAGAAGTTGCAGCTCGCCTTGACACAACAGGAAAAA
 CAAGAACAAGTTAAAATTTAGTCTCTAGCCAAGAGCTTAGAAGATGCTCTGAGGCAAACTGCAAGTGTCA
 CTCTGCAGGCTATTGCAGCTCAGAATGTGCGGTCCAGGCTGTCAATGCACACTCCAACATATTGAAGC
 CGCCATGGACAATTCTGAGATTGCAGGCGAGAAGAAATCTGCTCAGTGGCGCACAGTGGAGGGTGCATTG
 AAGGAACGCAGAAAGGCAGTAGATGAAGCTGCCGATGCCCTTCTCAAAGCCAAAGAAGAGTTAGAGAAGA
 TGAAAAGTGTGATTGAAAATGCAAAGAAAAAGAGGTTGCTGGGGCAAGCCTCATATACTGCTGCAGA
 GGGTAACTTTCACAACATGATAGTTGATCTGGATAATGTGGTCAAAAAGTCCAAGCAGCTCAGTCTGAG
 GCTAAGGTTGATCTCAGTATCATGAGCTGGTGGTCCAAGCTCGGGATGACTTAAACGAGAGCTGGACA
 GTATTACTCCAGAAGCTTCTGGGTGGAAAGGAATGAGTGTTCAGACTTAGCTGACAAGCTCTCTAC
 TGATGATCTGAACCTCCCTCATTGCTCATGCACATCGTCTGATTGATCAGCTGAACAGAGAGCTGGCAGAA
 CAGAAGGCCACCGAAAAGCAGCAGCATCACGTTAGCCTTGAGAAAACAAAAGCTGGAAGAAAAAGCGGCAT
 TTGACTCTGCAGTAGCAAAGCATTAGAACATCACAGAAGTGAATACAGGCTGAACAGGACAGAAAAGAT
 AGAAGAAGTCAGAGATGCCATGGAAAAATGAAATGAGAACCAGCTTCGCCGACAGGCAGCTGCCCACT
 GATCACTTGCAGATGTCCTTAGGGTACAAGAACAGGAATTGAAGTCTGAATTTGAGCAGAACCTGTCTG
 AGAACTCTCTGAACAAGAATTACAATTTGCTCGTCTCAGTCAAGAGCAAGTTGACAACCTTACTCTGGA
 TATAAATACTGCCTATGCCAGACTCAGAGGAATCGAACAGGCTGTTGAGGCCATGCAGTTGCTGAAGAG
 GAAGCCAGAAAAGCCACCAACTCTGGCTTTCAGTGGAGGCATTAAAGTACAGCATGAAGACCTCATCTG
 CAGAAACACCTACTATCCCGCTGGGTAGTGCAGTTGAGGCCATCAAAGCCAACTGTTCTGATAATGAATT
 CACCCAAGCTTTAACCGAGCTATCCCTCCAGAGTCCCTGACCCGTGGGGTGTACAGTGAAGAGACCCCT
 AGAGCCCGTTTCTATGCTGTTCAAAAATGGCCGAAAGGGTAGCAATGATTGATGAAACCAGAAATAGCT
 TGTACCACTACTTCTCTCTACCTACAGTCCCTGCTCCTATCCCACTCAGCAACTGAAGCCGCCCC
 AGAGCTCTGCCCTGAGGATATAAACACATTTAAATTAAGTGTATGCTTCTATTGATTGAGCATGGT
 GATCTGGAGCTAGCAGCAAAGTTTGTCAATCAGCTGAAGGGGGAATCCAGACGAGTGGCAGCAGGACTGGC
 TGAAGGAAGCCCGAATGACCCAGAAACGAAACAGATAGTGGAATCCTGACAGCATATGCCAGCGCCGT
 AGGAATAGGAACCACTCAGGTGCAGCCAGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC216114 representing NM_001100169
 Red=Cloning site Green=Tags(s)


MLRACQLSGVTAQAQSCLCGKFVLRPLRPCRRYSTSGSSGLTTGKIAGAGLLFVGGGIGGTILYAKWDSH
 FRESVEKTIPIYSDKLFEMVLGPAAYNVPLPKKSIQSGPLKISSVSEVMKESKQPASQLQKQKGDTPASAT
 APTEAAQIIISAAGDTLSVPAPAVQPEESLKTDPHEIGEGKPTALSEASSSIRERPPPEEVAARLAQOEK
 QEYQKIESLAKSLEDALRQTASVTLQAIQAQNAAVQAVNAHSNILKAAMDNSEIAGEKKSQAQWRTVEGAL
 KERRKAVDEAADALLKAKEELEKMKSVIENAKKKEVAGAKPHITAAEGKLNHMIVDLNDNVKQVAAQSE
 AKVVSQYHELVVQARDDFKRELDSTPEVLPGWKMSVSDLADKLSTDDLNSLIAHAHRRIDQLNRELAE
 QKATEKQHITLALQKLEEKRAFDSAVAKALEHHRSEIQAEQDRKIEEVRDAMENEMRTQLRRQAAAHT
 DHLRDVLRVQEQLKSEFEQNLSEKLEQELQFRRLSQEQVDFNLTLDINTAYARLRGIEQAVQSHAVAE
 EARKAHQLWLVEALKYSMTSSAETPTIPLGSAVEAIKANCSDNEFTQALTAAPPEL TRGVYSEETL
 RARFYAVQKLRVAMIDETRNLSYQYFLSYLQSLLLFPQQKPPPELCPEDINTFKLLSYASYCIEHG
 DLELAAKFVNQLKGESRRVAQDWLKEARMTLETQIVEILTAYASAVGIGTTQVQPE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001100169

ORF Size: 2271 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001100169.1](#), [NP_001093639.1](#)

RefSeq Size: 3015 bp

RefSeq ORF: 2274 bp

Locus ID: 10989


UniProt ID: [Q16891](#)

Cytogenetics: 2

MW: 83.5 kDa

Gene Summary: Component of the MICOS complex, a large protein complex of the mitochondrial inner membrane that plays crucial roles in the maintenance of crista junctions, inner membrane architecture, and formation of contact sites to the outer membrane. Plays an important role in the maintenance of the MICOS complex stability and the mitochondrial cristae morphology (PubMed:22114354, PubMed:25781180).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC216114