

Product datasheet for RC216155A1V

Human IL32 (NM_001012636) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human IL32 (NM_001012636) AAV Particle
Tag: Myc-DDK
Symbol: IL32
Synonyms: IL-32alpha; IL-32beta; IL-32delta; IL-32gamma; NK4; TAIF; TAIFa; TAIFb; TAIFc; TAIFd
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC216155 representing NM_001012636
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGGATCGCC

ATGTGCTTCCCGAAGGTCCTCTCTGATGACATGAAGAAGCTGAAGGCCCGAATGCACCAGGCCATAGAAA
 GATTTTATGATAAAATGCAAAATGCAGAATCAGGACGTGGACAGGACGACTTCAAAGAGGGCTACCTGGA
 GACAGTGGCGGCTTATTATGAGGAGCAGACCCAGAGCTCACTCCTCTACTTGAAAAAGAAAGAGATGGA
 TTACGGTGGCGAGGCAACAGATCCCCTGTCCCGGATGTTGAGGATCCCGCAACCGAGGAGCCTGGGGAGA
 GCTTTTGTGACAAGGTCATGAGATGGTCCAGGCCATGCTGCAGCGGCTGCAGACCTGGTGGCACGGGGT
 TCTGGCCTGGGTGAAGGAGAAGGTGGTGGCCCTGGTCCATGCAGTGCAGGCCCTCTGAAACAGTTCCAG
 AGTTTCTGCTGCTCTGTCAGAGCTTTCATGTCCTCTTTCCAGTCTACGGAGCCCCACGGGGGACA
 AGGAGGAGCTGACACCCAGAAGTGCTGAACCCCAATCCTCAAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC216155 representing NM_001012636
Red=Cloning site Green=Tags(s)

MCFPKVLSDDMKLKLKARMQAIERFYDKMQNAESGRGQDDFKEGYLETVAAYYEEQHPLETPLEKERDG
 LRCRGNRSPVPDVEDPATEEPGESFCDKVMRWFQAMLQRLQTWWHGVLAWVKEKVVAVHAVQALWKQFQ
 SFCCSLSELFMSSFQSYGAPRGDKEELTPQKCEPQSSK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



ACCN:	NM_001012636
ORF Size:	537 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001012636.1</u>
RefSeq Size:	880 bp
RefSeq ORF:	540 bp
Locus ID:	9235
UniProt ID:	<u>P24001</u>
Cytogenetics:	16p13.3
MW:	20.6 kDa