

Product datasheet for **RC216159L4V**

WSTF (BAZ1B) (NM_032408) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	WSTF (BAZ1B) (NM_032408) Human Tagged ORF Clone Lentiviral Particle
Symbol:	WSTF
Synonyms:	WBSCR9; WBSCR10; WSTF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032408
ORF Size:	4449 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216159).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032408.1 , NP_115784.1
RefSeq Size:	6043 bp
RefSeq ORF:	4452 bp
Locus ID:	9031
UniProt ID:	Q9UIG0
Cytogenetics:	7q11.23
Domains:	BROMO, PHD, DDT
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

MW: 170.7 kDa

Gene Summary: This gene encodes a member of the bromodomain protein family. The bromodomain is a structural motif characteristic of proteins involved in chromatin-dependent regulation of transcription. This gene is deleted in Williams-Beuren syndrome, a developmental disorder caused by deletion of multiple genes at 7q11.23. [provided by RefSeq, Jul 2008]