

Product datasheet for **RC216205L2V**

HSF4 (NM_001538) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | HSF4 (NM_001538) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | HSF4 |
| Synonyms: | CTM; CTRCT5 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_001538 |
| ORF Size: | 1386 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC216205). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001538.1 |
| RefSeq Size: | 1555 bp |
| RefSeq ORF: | 1389 bp |
| Locus ID: | 3299 |
| UniProt ID: | Q9ULV5 |
| Cytogenetics: | 16q22.1 |
| Domains: | HSF |
| Protein Families: | Druggable Genome, Transcription Factors |



[View online »](#)

MW: 49.8 kDa

Gene Summary: Heat-shock transcription factors (HSFs) activate heat-shock response genes under conditions of heat or other stresses. HSF4 lacks the carboxyl-terminal hydrophobic repeat which is shared among all vertebrate HSFs and has been suggested to be involved in the negative regulation of DNA binding activity. Two alternatively spliced transcripts encoding distinct isoforms and possessing different transcriptional activity have been described. [provided by RefSeq, Jul 2008]