

Product datasheet for **RC216218A1V**

Human WDR5 (NM_017588) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human WDR5 (NM_017588) AAV Particle
Tag:	Myc-DDK
Symbol:	WDR5
Synonyms:	BIG-3; CFAP89; SWD3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC216218 representing NM_017588 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGACGGAGGAGAAGAAGCCTGAGACCGAGGCCGCCAGAGCACAGCCAACCCCTTCGTCATCCGCCA
CTCAGAGCAAGCCTACACCTGTGAAGCCAACTATGCTCTAAAGTTCACCCTTGCTGGCCACACCAAAGC
AGTGTCCTCCGTGAAATTCAGCCCAATGGAGAGTGGCTGGCAAGTTCATCTGCTGATAAACTTATTTAA
ATTTGGGGCGCGTATGATGGGAAATTTGAGAAAACCATATCTGGTCACAAGCTGGGAATATCCGATGTAG
CCTGGTCGTGAGATTCTAACCTTCTTGTCTGCCTCAGATGACAAAACCTTGAAGATATGGGACGTGAG
CTCGGGCAAGTGTCTGAAAACCTGAAGGGACACAGTAATTATGTCTTTTGCTGCAACTTCAATCCCCAG
TCCAACCTTATTGTCTCAGGATCCTTTGACGAAAGCGTGAGGATATGGGATGTGAAAACAGGGAAGTGCC
TCAAGACTTTGCCAGCTCACTCGGATCCAGTCTCGGCCGTTCAATTTAATCGTGATGGATCCTTGATAGT
TTCAAGTAGCTATGATGGTCTCTGTGCGATCTGGGACACCGCCTCAGGCCAGTGCCTGAAGACGCTCATC
GATGACGACAACCCCCCGTGTCTTTGTGAAGTCTCCCCGAACGGCAAATACATCCTGGCCGCCACGC
TGGACAACACTCTGAAGCTCTGGGACTACAGCAAGGGGAAGTGCCTGAAGACGTACACTGGCCACAAGAA
TGAGAAATACTGCATATTTGCCAATTTCTGTACTGGTGGGAAGTGGATTGTGTCTGGCTCAGAGGAT
AACCTTGTGTTACATCTGGAACCTTCAGACGAAAGAGATTGTACAGAACTACAAGGCCACACAGATGTCTG
TGATCTCAACAGCTTGTACCCAAACAGAAAACATCATCGCCTCTGCTGCGCTAGAAAATGACAAAACAAT
TAAACTGTGGAAGAGTGACTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC216218 representing NM_017588 Red=Cloning site Green=Tags(s) MATEEKKPETEAARAQPTPSSSATQSKPTPVKPNYALKFTLAGHTKAVSSVKFSPNGEWLASSSADKLIK IWGAYDGKFEKTISGHKLGISDAWSSDSNLLVSASDDKTLKIWDVSSGKCLKTLKGHSNYVFCCNFPQ SNLIVSGSFDESRIWDVKTGKCLKTLPAHSDPVSAVHFNRDGLIVSSSYDGLCRIWDTASGQCLKTLI DDDNPPVSFVKFSPNGKYILAATLDNTLKLWDYSKGKCLKTYTGHKNEKCYCFANFSVTGGKWIVSGSED NLVYIWNLQTKETIVQKLQGHTDVIISTACHPTENIIASAALENDKTIKLWKSDC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_017588
ORF Size:	1002 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_017588.2</u>
RefSeq Size:	3163 bp
RefSeq ORF:	1005 bp
Locus ID:	11091
UniProt ID:	<u>P61964</u>
Cytogenetics:	9q34.2
MW:	36.4 kDa