

## Product datasheet for RC216235A1V

### Human ANKRD45 (NM\_198493) AAV Particle

#### Product data:

Product Type:	AAV Particles
Product Name:	Human ANKRD45 (NM_198493) AAV Particle
Tag:	Myc-DDK
Symbol:	ANKRD45
Synonyms:	CT117
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC216235 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**

ATGGAGTCAGAAGGACCTCCAGAGTCAGAGAGTTCAGAATTTTTCTCACAGCAAGAAGAGGAAAAATGAAGAAGAAGCCCAAGAACCAGAGGAAACAGGCCCTAAGAACCCTTTTACAACCTGCTCTCACAGGGGATGTAGAGGGTTTGCAGAAGATATTTGAGGATCCTGAGAATCCTCATCATGAACAGGCCATGCAGCTTCTTTAGAAGAAGACATCGTTGGGAGAAATTTGTTGTATGCAGCTTGCATGGCTGGGCAAAGTGACGTGATTGAGCTTTGGCAAATATGGTGTGAATCTGAATGAAAAACCACCAGAGGGTACACACTTTACATTTGCTGTCAGCCTGGGGTCGTTTTGGAACTTTGAAAGCACTGGTAGAACTGGATGTTGATATAGAAGCTTTGAACTTCCGGGAAGAAAGGGCTCGAGATGTTGCTGCTAGATATTCTCAGACTGAGTGTGTTGAATTCCTGGACTGGGCAGATGCAAGGCTGACTCTGAAAAATATATTGCAAAAGTCTCTTTAGCTGTTACTGACACAGAAAAAGGATCAGGAACTCCTTAAGGAAGACAAGAAATACCATCCTCAGTGCATGCCGTGCAAAAAATGAGTGGTTGAAACCCATACAGAAGCTTCCATTAATGAGCTTTTGGAGCAGAGACAACAACCTGGAAGATATTGTGATCCTATCTTCACAAAAATGACTACCCATGTCAAGTGAAGTGAAGTCCAAATCTGTAACAAGCCATGACCAAGAGAAGTCAAGATGATACCTCCAAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

<b>Protein Sequence:</b>	>RC216235 protein sequence Red=Cloning site Green=Tags(s)  MESEGPPESESSEFFSQQEENEEEEAEQEPEETGPKNPLLQPAL TGDVEGLQKIFEDPENPHHEQAMQLL LEEDIVGRNLLYAACMAGQSDVIRALAKYGVNLNEKTRGYTL LHCAAAGRLET LKALVELDVIDEALN FREERARDVAARYSQTECVFELDWADARL TLKKYIAKVSLAVTDTEKGSGLLKEDKNTILSACRAKNEW LETHTEASINELFEQRQLEDIVTPIFTKMTTPCQVKS AKSVTSHDQKRSQDDTSN  TRTRPLEQKLISEEDLAANDILDYKDDDDKV
<b>Species:</b>	Human
<b>Serotype:</b>	AAV-2
<b>ACCN:</b>	NM_198493
<b>ORF Size:</b>	798 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u>NM_198493.1</u>
<b>RefSeq Size:</b>	2671 bp
<b>RefSeq ORF:</b>	801 bp
<b>Locus ID:</b>	339416
<b>UniProt ID:</b>	<u>Q5TZF3</u>
<b>Cytogenetics:</b>	1q25.1
<b>MW:</b>	30 kDa