

## Product datasheet for **RC216246L4V**

### SLC6A15 (NM\_018057) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC6A15 (NM_018057) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC6A15
Synonyms:	hv7-3; NTT73; SBAT1; V7-3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018057
ORF Size:	867 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216246).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018057.4</a> , <a href="#">NP_060527.2</a>
RefSeq Size:	4227 bp
RefSeq ORF:	870 bp
Locus ID:	55117
UniProt ID:	<a href="#">Q9H2J7</a>
Cytogenetics:	12q21.31
Domains:	SNF
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

**MW:** 32.1 kDa

**Gene Summary:** This gene encodes a member of the solute carrier family 6 protein family which transports neutral amino acids. The encoded protein is thought to play a role in neuronal amino acid transport (PMID: 16185194) and may be associated with major depression (PMID: 21521612). Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2012]