

Product datasheet for **RC216270L3V**

REXO1L1 (REXO1L1P) (NM_172239) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	REXO1L1 (REXO1L1P) (NM_172239) Human Tagged ORF Clone Lentiviral Particle
Symbol:	REXO1L1
Synonyms:	GOR; REXO1L1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_172239
ORF Size:	2025 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216270).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172239.2 , NP_758439.4
RefSeq Size:	7032 bp
RefSeq ORF:	2027 bp
Locus ID:	254958
Cytogenetics:	8q21.2
MW:	73.7 kDa



[View online »](#)