

Product datasheet for **RC216294L2V**

ARID5A (NM_212481) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ARID5A (NM_212481) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ARID5A
Synonyms:	MRF-1; MRF1; RP11-363D14
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_212481
ORF Size:	1782 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216294).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_212481.1
RefSeq Size:	2223 bp
RefSeq ORF:	1785 bp
Locus ID:	10865
UniProt ID:	Q03989
Cytogenetics:	2q11.2
MW:	63.9 kDa



[View online »](#)

Gene Summary:

Members of the ARID protein family, including ARID5A, have diverse functions but all appear to play important roles in development, tissue-specific gene expression, and regulation of cell growth (Patsialou et al., 2005 [PubMed 15640446]).[supplied by OMIM, Mar 2008]