

Product datasheet for **RC216303A1V**

Human **RAB6B** (NM_016577) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RAB6B (NM_016577) AAV Particle
Tag: Myc-DDK
Symbol: RAB6B
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC216303 representing NM_016577
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCGCGAGGGGAGATTTTGGGAATCCACTGAGAAAATTCAAGTTGGTGTCTTGGGGAGCAGAGCG
TCGGGAAGACGTCTCTGATTACGAGGTTTCATGTACGACAGCTTCGACAACACATACCAGGCAACCATTGG
GATTGACTTCTGTCAAAAACCATGTACTTGGAGGACCGCACGGTCCGACTGCAGCTCTGGGACACAGCT
GGTCAGGAGAGGTTCCGAGCCTGATCCCCAGCTACATCCGGGACTCCACGGTGGCTGTGGTGGTGTACG
ACATCACAAATCTCAACTCCTTCCAACAGACCTCTAAGTGGATCGACGACGTGAGGACAGAGAGGGGCAG
TGATGTTATCATCATGCTGGTGGCAACAAGACGGACCTGGCTGATAAGAGGCAGATAAACCATCGAGGAG
GGGAGCAGCGCGCAAGAAGAACTGAGCGTCATGTTTATTGAGACCAGTGCAGAACTGGCTACAACGTGA
AGCAGCTTTTTCGACGTGTGGCGTCGGCTCTACCCGGAATGGAGAATGTCCAGGAGAAAAGCAAAGAAGG
GATGATTGACATCAAGCTGGACAAACCCAGGAGCCCCCGGCCAGCGAGGGCGGCTGCTCCTGC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC216303 representing NM_016577
Red=Cloning site Green=Tags(s)

MSAGGDFGNPLRKFLVFLGEQSVGKTSITRFMYDSFDNTYQATIGIDFLSKTMYLEDRTVRLQLWDTA
GQERFRSLIPSYIRDSTVAVVVYDITNLNSFQQTSKWIDDVRTERGSDVIIMLVGNKTDLADKRQITIEE
GEQRAKELSVMFIE TSAKTGYNVKQLFRRVASALPGMENVQESKEGMIDIKLDKPQEPPESEGGCSC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_016577
ORF Size:	624 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_016577.3</u>
RefSeq Size:	5561 bp
RefSeq ORF:	627 bp
Locus ID:	51560
UniProt ID:	<u>Q9NRW1</u>
Cytogenetics:	3q22.1
MW:	23.3 kDa