

Product datasheet for **RC216330L3V**

SLFN5 (NM_144975) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLFN5 (NM_144975) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLFN5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144975
ORF Size:	2673 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216330).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144975.3
RefSeq Size:	4660 bp
RefSeq ORF:	2676 bp
Locus ID:	162394
UniProt ID:	Q08AF3
Cytogenetics:	17q12
Domains:	AAA_div
MW:	101 kDa
Gene Summary:	May have a role in hematopoietic cell differentiation.[UniProtKB/Swiss-Prot Function]



[View online »](#)