

Product datasheet for **RC216339L3V**

AVL9 (NM_015060) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AVL9 (NM_015060) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AVL9
Synonyms:	KIAA0241
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_015060
ORF Size:	1944 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216339).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015060.1 , NP_055875.1
RefSeq Size:	2353 bp
RefSeq ORF:	1947 bp
Locus ID:	23080
UniProt ID:	Q8NBF6
Cytogenetics:	7p14.3
MW:	71.8 kDa
Gene Summary:	Functions in cell migration.[UniProtKB/Swiss-Prot Function]


[View online »](#)