

Product datasheet for RC216355A1V

Human RPL30 (NM_000989) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RPL30 (NM_000989) AAV Particle
Tag: Myc-DDK
Symbol: RPL30
Synonyms: L30
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC216355 representing NM_000989
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTGGCCGCAAGAAGACGAAAAAGTCGCTGGAGTCGATCAACTCTAGGCTCCAACCTCGTTATGAAAA
 GTGGGAAGTACGTCCTGGGGTACAAGCAGACTCTGAAGATGATCAGACAAGGCAAAGCGAAATTTGGTCAT
 TCTCGTAACAACCTGCCAGCTTTGAGGAAATCTGAAATAGAGTACTATGCTATGTTGGCTAAAACCTGGT
 GTCCATCACTACAGTGGCAATAATTTGAACTGGGCACAGCATGCGGAAAATACTACAGAGTGTGCACAC
 TGGCTATCATTGATCCAGGTGACTCTGACATCATTAGAAGCATGCCAGAACAGACTGGTAAAAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC216355 representing NM_000989
 Red=Cloning site Green=Tags(s)

MVAAKKTKKSLESINSRLQLVMKSGKYVLGYKQTLKMIROGKAKLVILANNCPALRKSEIEYYAMLAKTG
 VHHYSGNNIELGTACGKYRVCTLAIIDPGSDIIRSMPEQTGEK

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_000989
ORF Size: 345 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_000989.2</u>
RefSeq Size:	524 bp
RefSeq ORF:	348 bp
Locus ID:	6156
UniProt ID:	<u>P62888</u>
Cytogenetics:	8q22.2
MW:	12.6 kDa