

## Product datasheet for **RC216366L4V**

### PPFIA4 (NM\_015053) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PPFIA4 (NM_015053) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PPFIA4
Synonyms:	Liprin-alpha4; liprin alpha4; protein tyrosine phosphatase, receptor type, f polypeptide (PTPRF), interacting protein (liprin), alpha 4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015053
ORF Size:	2103 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216366).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_015053.1</a> , <a href="#">NP_055868.1</a>
RefSeq Size:	5055 bp
RefSeq ORF:	2106 bp
Locus ID:	8497
Cytogenetics:	1q32.1
Protein Families:	Druggable Genome, Phosphatase
MW:	77.9 kDa


[View online »](#)

**Gene Summary:**

PPFIA4, or liprin-alpha-4, belongs to the liprin-alpha gene family. See liprin-alpha-1 (LIP1, or PPFIA1; MIM 611054) for background on liprins.[supplied by OMIM, Mar 2008]