

Product datasheet for **RC216382L4V**

ZNF286 (ZNF286A) (NM_020652) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF286 (ZNF286A) (NM_020652) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF286
Synonyms:	ZNF286
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020652
ORF Size:	1563 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216382).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020652.1
RefSeq Size:	2515 bp
RefSeq ORF:	1566 bp
Locus ID:	57335
UniProt ID:	Q9HBT8
Cytogenetics:	17p12
Domains:	KRAB, zf-C2H2
Protein Families:	Transcription Factors



[View online »](#)

MW: 60 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]