

Product datasheet for **RC216399L3V**

ERC1 (NM_015064) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ERC1 (NM_015064) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ERC1
Synonyms:	Cast2; ELKS; ELKS/RAB6-interacting/CAST family member 1; KIAA1081; MGC12974; OTTHUMP00000160633; RAB6 interacting protein 2; RAB6-interacting protein 2; RAB6IP2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_015064
ORF Size:	2844 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216399).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015064.2 , NP_055879.1
RefSeq Size:	6919 bp
RefSeq ORF:	2846 bp
Locus ID:	23085
Cytogenetics:	12p13.33
Domains:	Pox_A_type_inc
Protein Families:	Druggable Genome



[View online »](#)

MW: 108.6 kDa

Gene Summary: The protein encoded by this gene is a member of a family of RIM-binding proteins. RIMs are active zone proteins that regulate neurotransmitter release. This gene has been found fused to the receptor-type tyrosine kinase gene RET by gene rearrangement due to the translocation t(10;12)(q11;p13) in thyroid papillary carcinoma. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]