

Product datasheet for **RC216415L4V**

RTN4RL1 (NM_178568) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RTN4RL1 (NM_178568) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RTN4RL1
Synonyms:	NgR3; NGRH2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_178568
ORF Size:	1323 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216415).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178568.2
RefSeq Size:	3183 bp
RefSeq ORF:	1326 bp
Locus ID:	146760
UniProt ID:	Q86UN2
Cytogenetics:	17p13.3
Protein Families:	Druggable Genome
MW:	49.1 kDa



[View online »](#)

Gene Summary:

Cell surface receptor. Plays a functionally redundant role in postnatal brain development and in regulating axon regeneration in the adult central nervous system. Contributes to normal axon migration across the brain midline and normal formation of the corpus callosum. Protects motoneurons against apoptosis; protection against apoptosis is probably mediated by MAG. Plays a role in inhibiting neurite outgrowth and axon regeneration via its binding to neuronal chondroitin sulfate proteoglycans. Binds heparin (By similarity). Like other family members, plays a role in restricting the number dendritic spines and the number of synapses that are formed during brain development (PubMed:22325200). Signaling mediates activation of Rho and downstream reorganization of the actin cytoskeleton (PubMed:22325200). [UniProtKB/Swiss-Prot Function]