

Product datasheet for **RC216428**

FGFR3 (NM_022965) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FGFR3 (NM_022965) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	FGFR3
Synonyms:	ACH; CD333; CEK2; HSGFR3EX; JTK4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC216428 representing NM_022965
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGCGCCCCTGCCTGCGCCCTCGCGCTCTCGTGGCCGTGGCCATCGTGGCCGGCGCCTCCTCGGAGT
 CCTTGGGACGGAGCAGCGCGTCTGGGGCAGCGGCAGAAGTCCCGGGCCAGAGCCCGCCAGCAGGA
 GCAGTTGGTCTTCGGCAGCGGGATGCTGTGGAGCTGAGCTGTCCCCGCCCCGGGGTGGTCCCATGGGG
 CCCACTGTCTGGGTCAAGGATGGCACAGGGCTGGTGCCCTCGGAGCGTGTCTGGTGGGGCCCCAGCGGC
 TGCAGGTGCTGAATGCCTCCACGAGGACTCCGGGGCCTACAGCTGCCGGCAGCGGCTCACGCAGCGCT
 ACTGTGCCACTTCAGTGTGCGGGTGACAGACGCTCCATCCTCGGGAGATGACGAAGACGGGGAGGACGAG
 GCTGAGGACACAGGTGTGGACACAGGGGCCCTTACTGGACACGGCCGAGCGGATGGACAAGAAGCTGC
 TGGCCGTGCCGGCCCAACACCGTCCGCTTCCGCTGCCAGCCGCTGGCAACCCCACTCCCTCCATCTC
 CTGGCTGAAGAACGGCAGGGAGTTCGCGGCGAGCACCGCATTGGAGGCATCAAGCTGCGGCATCAGCAG
 TGGAGCCTGGTCATGAAAGCGTGGTGCCCTCGGACCGCGCAACTACACCTGCGTCTGGAGAACAAGT
 TTGGCAGCATCCGGCAGACGTACACGCTGGACGTGCTGGAGCGCTCCCCGACCGGCCATCCTGCAGGC
 GGGGCTGCCGGCAACCAGACGGCGGTGCTGGGCAGCGACGTGGAGTTCCACTGCAAGGTGTACAGTGAC
 GCACAGCCCCACATCCAGTGGCTCAAGCAGTGGAGGTGAACGGCAGCAAGGTGGGCCGGACGGCACAC
 CCTACGTTACCGTGTCAAGGTGTCCTGGAGTCCAACCGTCCATGAGCTCCAACACACCACTGGTGCG
 CATCGAAGGCTGTCTCAGGGGAGGGCCCCACGCTGGCCAATGTCTCCGAGCTCGAGTGCCTGCCGAC
 CCCAAATGGGAGCTGTCTCGGGCCCGCTGACCCTGGCAAGCCCTTGGGGAGGGCTGCTTCGGCCAGG
 TGGTCATGGCGGAGGCCATCGGCATTGACAAGGACCGGGCCCAAGCCTGTACCGTAGCCGTGAAGAT
 GCTGAAAGACGATGCCACTGACAAGGACCTGTCGGACCTGGTGTCTGAGATGGAGATGATGAAGATGATC
 GGGAAACACAAAAACATCATCAACCTGCTGGGCGCCTGCACGCAGGGCGGGCCCCGTACGTGCTGGTGG
 AGTACGCGGCCAAGGTAACCTGCGGGAGTTTCTGCGGGCGCGCGGCCCGGGCTGGACTACTCCTT
 CGACACCTGCAAGCCGCCGAGGAGCAGCTCACCTTCAAGGACCTGGTGTCTGTGCCTACCAGGTGGCC
 CGGGGCATGGAGTACTTGGCCTCCAGAAGTGCATCCACAGGGACCTGGCTGCCCGCAATGTGCTGGTGA
 CCGAGGACAACGTGATGAAGATCGCAGACTTCGGGCTGGCCGGGACGTGCACAACCTCGACTACTACAA
 GAAGACAACCAACGGCCGGCTGCCCGTGAAGTGGATGGCGCTGAGGCCCTGTTTGACCGAGTCTACACT
 CACCAGAGTGACGTCTGGTCTTTGGGGTCTGCTCTGGGAGATCTTACGCTGGGGGGCTCCCCGTACC
 CCGGCATCCCTGTGGAGGAGCTTCAAGCTGCTGAAGGAGGGCCACCGCATGGACAAGCCCGCCAACCTG
 CACACACGACCTGTACATGATCATGCGGGAGTGTGGCATGCCGCGCCCTCCAGAGGCCACCTTCAAG
 CAGCTGGTGGAGGACCTGGACCGTGTCTTACCCTGACGTCCACCGACGAGTACCTGGACCTGTCCGGCGC
 CTTTCGAGCAGTACTCCCGGGTGGCCAGGACACCCCAAGCTCCAGCTCCTCAGGGGACGACTCCGTGTT
 TGCCACGACCTGCTGCCCGGGCCCCACCCAGCAGTGGGGGCTCGCGGACG

ACGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGA
TTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC216428 representing NM_022965
Red=Cloning site Green=Tags(s)


```

MGAPACALALCVAVAIVAGASSESLGTEQRVVGRAAEVPPEPGQEQQLVFGSGDAVELSCPPPGGGPMG
PTVWVKDGTGLVPSERVLVGPQRLQVLNASHEDSGAYSCRQRLTQRVLCFVSVRVTDAPSSGDDDEGEDE
AEDTGVDTGAPYWTRPERMDKLLAVPAANTVRFRCPAAGNTPPSISWLKNGREFRGEHRIGGIKLRHQQ
WSLVMESVVPSSDRGNYTCVVENKFGSIRQTYTLDVLESPHRPILQAGLPANQTAVLGSDVEFHCKVYSD
AQPHIQWLKHVEVNGSKVGPDPYVTVLKVLSLESNASMSNTPLVRIARLSSGEGPTLANVSELELPAD
PKWELSRARLTLGKPLGEGCFGQVYMAEAIIGIDKDRAAKPVTAVKMLKDDATDKDLSLVSEMEMMKMI
GKHKNIINLLGACTQGGPLYVLYEYAAKGNLREFLRARRPPGLDYSDFTCKPPEEQQLTFKDLVSCAYQVA
RGMEYLASQKCIHRDLAARNVLVTDENVMKIADFLARDVHNLDYKKTNGRLPVKWMPEALFDRVYT
HQSDVWSFGVLLWEIFTLGGSPYGPVPEELFKLLKEGHRMDKPANCTHDLYMIMRECWHAAPSQRPTFK
QLVEDLDRVLTVTSTDEYLDLAPFEEQYSPGGQTPSSSSSGDDSVFAHDLLPPAPPSSGGSR
    
```

TRRLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-NotI

Cloning Scheme:



ACCN: NM_022965

ORF Size: 2082 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_022965.1 , NP_075254.1
RefSeq Size:	3757 bp
RefSeq ORF:	2085 bp
Locus ID:	2261
UniProt ID:	P22607
Cytogenetics:	4p16.3
Domains:	pkinase, TyrKc, S_TKc, ig, IGc2, IG
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane
Protein Pathways:	Bladder cancer, Endocytosis, MAPK signaling pathway, Pathways in cancer, Regulation of actin cytoskeleton
MW:	73.7 kDa
Gene Summary:	<p>This gene encodes a member of the fibroblast growth factor receptor (FGFR) family, with its amino acid sequence being highly conserved between members and among divergent species. FGFR family members differ from one another in their ligand affinities and tissue distribution. A full-length representative protein would consist of an extracellular region, composed of three immunoglobulin-like domains, a single hydrophobic membrane-spanning segment and a cytoplasmic tyrosine kinase domain. The extracellular portion of the protein interacts with fibroblast growth factors, setting in motion a cascade of downstream signals, ultimately influencing mitogenesis and differentiation. This particular family member binds acidic and basic fibroblast growth hormone and plays a role in bone development and maintenance. Mutations in this gene lead to craniosynostosis and multiple types of skeletal dysplasia. [provided by RefSeq, Aug 2017]</p>

Product images:



Circular map for RC216428