

Product datasheet for **RC216463L2V**

MAP3K12 (NM_006301) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MAP3K12 (NM_006301) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MAP3K12
Synonyms:	DLK; MEKK12; MUK; ZPK; ZPKP1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_006301
ORF Size:	2577 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216463).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006301.3
RefSeq Size:	3541 bp
RefSeq ORF:	2580 bp
Locus ID:	7786
UniProt ID:	Q12852
Cytogenetics:	12q13.13
Domains:	pkinese, TyrKc, S_TKc
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

Protein Pathways: MAPK signaling pathway

MW: 93.7 kDa

Gene Summary: This gene encodes a member of the serine/threonine protein kinase family. This kinase contains a leucine-zipper domain and is predominately expressed in neuronal cells. The phosphorylation state of this kinase in synaptic terminals was shown to be regulated by membrane depolarization via calcineurin. This kinase forms heterodimers with leucine zipper containing transcription factors, such as cAMP responsive element binding protein (CREB) and MYC, and thus may play a regulatory role in PKA or retinoic acid induced neuronal differentiation. Alternatively spliced transcript variants encoding different proteins have been described.[provided by RefSeq, Jul 2010]