

Product datasheet for **RC216473L4V**

MAEA (NM_005882) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MAEA (NM_005882) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MAEA
Synonyms:	EMLP; EMP; GID9; HLC-10; P44EMLP; PIG5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005882
ORF Size:	1065 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216473).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005882.2
RefSeq Size:	1983 bp
RefSeq ORF:	1068 bp
Locus ID:	10296
UniProt ID:	Q7L5Y9
Cytogenetics:	4p16.3
Domains:	LisH, CTLH
Protein Families:	Druggable Genome



[View online »](#)

MW: 40.2 kDa

Gene Summary: This gene encodes a protein that mediates the attachment of erythroblasts to macrophages. This attachment promotes terminal maturation and enucleation of erythroblasts, presumably by suppressing apoptosis. The encoded protein is an integral membrane protein with the N-terminus on the extracellular side and the C-terminus on the cytoplasmic side of the cell. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]