

## Product datasheet for **RC216493**

### OLFML2A (NM\_182487) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	OLFML2A (NM_182487) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	OLFML2A
Synonyms:	PRO34319
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RC216493 representing NM\_182487  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCCGCTGCCGCCCTCCCGCCCGGCCGCTGCTCCTTCTGCCGCTAGTGCTGCTGCTGAGCGGCCGCC  
 CCACGCGCGCCGACAGTAAGGTGTTGGGGACCTGGACCAGTGAGGATGACCTCGGAGGGCTCCGACTG  
 CCGTTGCAAGTGCATCATGCGGCCCTGAGCAAGGACGCGTGTAGCCGAGTGCAGTGGCGGGCAGCG  
 GTGGAGGACTTCTACACGGTGGAGACTGTGAGCTCGGGCACTGACTGCCGCTGCTCCTGTACCGCACCTC  
 CCTCCTCTCTCAACCCCTGTGAGAACGAGTGGAAGATGGAGAACTCAAAAAGCAGGCGCCCGAGCTCCT  
 CAAGCTGCAGTCCATGGTGGATCTCCTGGAGGGCACCTGTACAGCATGGACTTGTGAAGTGCACGCC  
 TACGTCCACAAGGTGGCCTCCAGATGAACACACTGGAAGAGAGCATCAAGGCCAACCTGAGCCGGGAGA  
 ATGAGGTGGTGAAGGACAGCGTGCACCCTCAGTGAGCAGTTGAGGCACTATGAGAATCACTCTGCCAT  
 CATGCTGGGCATCAAGAAGGAGCTGTCCCGCTGGGCCTCCAGCTGCTGCAGAAGGATGCCGCCCGGCC  
 CCTGCCACCCCTGCCACGGGCACTGGTAGCAAGGCCAGGACACAGCTAGAGGAAAAGGCCAAGGACATCA  
 GCAAGTATGGCAGTGTGCAGAAAAGCTTTCAGACAGAGGCCCTCCAAAACCTCCCAAGGAGAAGCTGCT  
 TCAGGTGGAGAAGCTGAGAAAAGGAGAGCGGCAAGGGCAGTTTCTCCAGCCACAGCCAAGCCCCGCGCC  
 CTGGCCAGCAGCAGGCTGTGATCCGGGGCTTACCTACTACAAGGCAGGCAAGCAGGAGGTGACCGAGG  
 CGGTGGCAGACAACCCCTCCAGGGCACTTCTGGCTGGAGCAACTGCCGCCAAGGTGGAGGGCAGGTC  
 CAACTCCGAGAGCCCAACTCCGAGAGCAGGATGAGGCTGAGCCAGGTCTCCGAGCGAGTGGACCTG  
 GCTTCTGGCACCCCACTTCAATCCCTGCCACCACCACCACCACCACCACCCCAACCCCAACCCCAACCA  
 GTCTCTGCCACCGAGCCACTTCAGGTCCAGAAGTCTCCAGCCAAGGCAGAGAGGCGAGCTGTGAGGG  
 CACCCTCCGGGCTGTGGACCCCTGTGAGGCACCACAGCTATGGGCGCCACGAGGGAGCCTGGATGAAG  
 GACCCTGCAGCTCGAGACGACAGGATCTATGTCACCAACTACTACTATGGAAACAGCCTGGTGGAGTTCC  
 GCAACCTGAAAACTTCAAGCAAGGCCGCTGGAGTAACATGTACAAGCTACCCTACAAGTGGATCGGCAC  
 AGGCCACGTGGTGTACCAGGGCGCCTTCTACTACAACCGCGCCTTACCAAGAACATCATCAAGTACGAC  
 CTACGGCAGCGCTTCGTGGCCTCCTGGGCGCTGCTGCCGACGTGGTATATGAGGACACCACACCTTGA  
 AGTGGCGCGGACACTCGGACATTGACTTTGCCGTGGACGAGAGCGGCCTGTGGGTCTCTACCCCGCGT  
 GGACGACCGGATGAGGCCAGCCGAGGTGATCGTCTGAGTCGCTTGGACCCCGGCGATCTCTCCGTG  
 CACCGGGAGACCAGTGAAGACACGGCTGCGGCGGAACCTCTACGGGAACCTGCTTCTGGTGTGCGGCA  
 TCCTGTATGCCGTGGACACGTACAACCAGCAGGAAGGCCAGGTGCGCTACGCTTTCGACACGCACACGGG  
 CACCGACGCACGCCCCAGCTGCCGTTCTCAACGAGCAGCCTACACCACCCAGATCGACTACAACCC  
 AAGGAGCGGGTGTGTACGCCTGGGACAATGGCCACCAGCTCACCTACCCCTCCACTTCGTGGTC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGATAAGGTTAA

**Protein Sequence:** >RC216493 representing NM\_182487  
Red=Cloning site Green=Tags(s)

```


MAAAALPPRPLLLLPLVLLLSGRPTRADSKVFGDLQVRMTSEGSDCRCKCIMRPLSKDACSRVRSGRAR
VEDFYTVETVSSGTDRCRCSTAPPSSLNPCENEWKMEKLLKQAPPELLKLSMVDLLEGLYSMDLMKVHA
YVHKVASQMNTLEESIKANLSRENEVVKDSVRHLSEQLRHYENHSAIMLGIKKELSRLLGLQLLQKDAAAA
PATPATGTGSKAQDTARGKGGKDISKYGSVQKSFADRLPKPPKEKLLQVEKLRKESGKGSFLQPTAKPRA
LAQQQAVIRGFYYKAGKQEVTEAVADNTLQGTSWLEQLPPKVEGRSNSAEPNSAEQDEAEPRSSSERVDL
ASGTPTSIPATTTTATTTPTPTTSLPTEPPSGPEVSSQGREASCEGTLRAVDPPVRHHSYGRHEGAWMK
DPAARDDRIYVTNYYYGNLVEFRNLENFKQGRWSNMYKLPYNWIGTGHVYVYQGAFFYNRAFTKNIKYD
LRQRFVSWALLPDVVYEDTTPWKWRGHSIDFAVDESLWVIYPAVDDRDEAQPEVIVLSRLDPGDL SV
HRETTWKTRLRNSYGNCFVCGILYAVDTYNQEQGVAYAFDTHGTDARPLPFLNEHAYTTQIDYNP
KERVLYAWDNGHQLTYTLHFVV
    
```

SGPTRRRLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk8016\\_g09.zip](https://cdn.origene.com/chromatograms/mk8016_g09.zip)

**Restriction Sites:** SgfI-RsrII

**Cloning Scheme:**



**ACCN:** NM\_182487

**ORF Size:** 1956 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_182487.4](#)

**RefSeq Size:** 6556 bp

**RefSeq ORF:** 1959 bp

**Locus ID:** 169611


**UniProt ID:** [Q68BL7](#)

**Cytogenetics:** 9q33.3

**Protein Families:** Druggable Genome

**MW:** 72.9 kDa

**Product images:**



Circular map for RC216493