

Product datasheet for **RC216508L3V**

CNG3 (CNGA3) (NM_001079878) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CNG3 (CNGA3) (NM_001079878) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CNG3
Synonyms:	ACHM2; CCNC1; CCNCa; CCNCalpha; CNGG3; CNG3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001079878
ORF Size:	2028 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216508).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001079878.1 , NP_001073347.1
RefSeq Size:	3794 bp
RefSeq ORF:	2031 bp
Locus ID:	1261
UniProt ID:	Q16281
Cytogenetics:	2q11.2
Protein Families:	Druggable Genome, Ion Channels: Cyclic nucleotide gated, Transmembrane
Protein Pathways:	Olfactory transduction



[View online »](#)

MW: 76.9 kDa

Gene Summary: This gene encodes a member of the cyclic nucleotide-gated cation channel protein family which is required for normal vision and olfactory signal transduction. Mutations in this gene are associated with achromatopsia (rod monochromacy) and color blindness. Two alternatively spliced transcripts encoding different isoforms have been described. [provided by RefSeq, Jul 2008]