

## Product datasheet for **RC216534L3V**

### HPSE2 (NM\_021828) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HPSE2 (NM_021828) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HPSE2
Synonyms:	HPA2; HPR2; UFS; UFS1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_021828
ORF Size:	1776 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216534).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021828.3</a>
RefSeq Size:	4349 bp
RefSeq ORF:	1779 bp
Locus ID:	60495
UniProt ID:	<a href="#">Q8WWQ2</a>
Cytogenetics:	10q24.2
Protein Pathways:	Glycosaminoglycan degradation, Metabolic pathways
MW:	66.6 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a heparanase enzyme. The encoded protein is an endoglycosidase that degrades heparin sulfate proteoglycans located on the extracellular matrix and cell surface. This protein may be involved in biological processes involving remodeling of the extracellular matrix including angiogenesis and tumor progression. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]