

Product datasheet for RC216590L2V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CD45 (PTPRC) (NM_002838) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: CD45 (PTPRC) (NM 002838) Human Tagged ORF Clone Lentiviral Particle

Symbol: CD45

Synonyms: B220; CD45; CD45R; GP180; L-CA; LCA; LY5; T200

Mammalian Cell

Selection:

None

Vector: pLenti-C-mGFP (PS100071)

Tag: mGFP

ACCN: NM_002838 **ORF Size:** 3918 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC216590).

Sequence:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 002838.3

 RefSeq Size:
 5026 bp

 RefSeq ORF:
 3921 bp

 Locus ID:
 5788

 UniProt ID:
 P08575

Cytogenetics: 1q31.3-q32.1

Domains: Y_phosphatase, PTPc_motif, FN3

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Phosphatase, Transmembrane





CD45 (PTPRC) (NM_002838) Human Tagged ORF Clone Lentiviral Particle - RC216590L2V

Protein Pathways: Cell adhesion molecules (CAMs), Fc gamma R-mediated phagocytosis, Primary

immunodeficiency, T cell receptor signaling pathway

MW: 147.3 kDa

Gene Summary: The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP)

family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitosis, and oncogenic transformation. This PTP contains an extracellular domain, a single transmembrane segment and two tandem

intracytoplasmic catalytic domains, and thus is classified as a receptor type PTP. This PTP has

been shown to be an essential regulator of T- and B-cell antigen receptor signaling. It functions through either direct interaction with components of the antigen receptor complexes, or by activating various Src family kinases required for the antigen receptor signaling. This PTP also suppresses JAK kinases, and thus functions as a regulator of cytokine receptor signaling. Alternatively spliced transcripts variants of this gene, which encode

distinct isoforms, have been reported. [provided by RefSeq, Jun 2012]