

Product datasheet for **RC216610L3V**

TFAP2D (NM_172238) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | TFAP2D (NM_172238) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | TFAP2D |
| Synonyms: | TFAP2BL1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_172238 |
| ORF Size: | 1356 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC216610). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_172238.1 |
| RefSeq Size: | 2040 bp |
| RefSeq ORF: | 1359 bp |
| Locus ID: | 83741 |
| UniProt ID: | Q7Z6R9 |
| Cytogenetics: | 6p12.3 |
| Protein Families: | Transcription Factors |
| MW: | 49.4 kDa |



[View online »](#)

Gene Summary:

Sequence-specific DNA-binding protein that interacts with inducible viral and cellular enhancer elements to regulate transcription of selected genes. AP-2 factors bind to the consensus sequence 5'-GCCNNNGGC-3' and activate genes involved in a large spectrum of important biological functions including proper eye, face, body wall, limb and neural tube development. They also suppress a number of genes including MCAM/MUC18, C/EBP alpha and MYC (By similarity).[UniProtKB/Swiss-Prot Function]