

Product datasheet for **RC216623A1V**

Human Folate Receptor 4 (IZUMO1R) (NM_001080486) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Folate Receptor 4 (IZUMO1R) (NM_001080486) AAV Particle
Tag:	Myc-DDK
Symbol:	Folate Receptor 4
Synonyms:	folate receptor 4 (delta) homolog; folate receptor 4 (delta) homolog (mouse); Folbp3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC216623 representing NM_001080486 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCATGCTGGTGGCCGCTCCTGCTAGAGCTGTGGACAGTCATGCCACCTGGGCTGGGGACGAGCTGC
TCAACATCTGCATGAATGCCAAACACCACAAGAGAGTGCCAGCCAGCCAGAAGACAAGCTCTATGAGGAGTG
CATCCCCTGGAAGGACAATGCCTGTGCACCCTCACGACAAGCTGGGAAGCCCATCTGGATGTATCCCCA
CTCTACAACTTCAGCCTGTTTCACTGTGGACTGCTGATGCCTGGCTGTCGGAAGCACTTCATCCAGGCTA
TCTGCTTCTATGAGTGCTCCCAAACCTGGGGCCCTGGATCCAGCCAGTGCCCGGAGTGCGCAGGGAGA
GCGAGTTGTGAATGTGCCGCTGTGCCAGGAGGACTGTGAGGAGTGGTGGGAAGACTGTCCGATGTCTTAC
ACATGCAAAATCCAATGGCGTGGTGGCTGGGACTGGAGTCAGGGGAAGAACCGCTGCCCCAAAGGGGCC
AGTGCCTCCCTTTCTCCATTACTTCCCCACCCAGCTGACCTGTGTGAGAAGACTTGGAGCAATTCCTT
CAAAGCCAGCCCTGAGCGACGGAACAGTGGCGGTGTCTCCAGAAGTGGTTTGAGCCTGCTCAGGCAAC
CCCAATGTGGCCGTGGCCGCTCTTCGCCAGCTCTGCCCCATCTGGGAAGTGTCTACACCATCATGG
TCTGCTCCCTGTTCTGCCGTTCTTTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC216623 representing NM_001080486 Red=Cloning site Green=Tags(s) MACWWPLLELWTVMPTWAGDELLNICMNAKHHKRVSPEDKLYEECIPWKDNACCTLTTSWEAHLDVSP LYNFSLFHCGLLMPGCRKHFIQAICFYECSPNLGPWIQPVAPSGQGERVVNVPLCQEDCEEWWEDCRMSY TCKSNWRGGWDWSQGNRCPKGAQCLPFSHYFPTPADLCEKTWSNSFKASPERRNSGRCLQKWFEP AQGN PNVAVARLFASSAPSWELSYTIMVCSLFLPFLS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001080486
ORF Size:	729 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001080486.1, NP_001073955.1</u>
RefSeq Size:	732 bp
RefSeq ORF:	731 bp
Locus ID:	390243
Cytogenetics:	11q21
MW:	27.8 kDa