

Product datasheet for **RC216626L3V**

BICC1 (NM_001080512) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BICC1 (NM_001080512) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BICC1
Synonyms:	BICC; CYSRD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001080512
ORF Size:	2922 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216626).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001080512.1
RefSeq Size:	3119 bp
RefSeq ORF:	2925 bp
Locus ID:	80114
UniProt ID:	Q9H694
Cytogenetics:	10q21.1
MW:	104.7 kDa



[View online »](#)

Gene Summary:

This gene encodes an RNA-binding protein that is active in regulating gene expression by modulating protein translation during embryonic development. Mouse studies identified the corresponding protein to be under strict control during cell differentiation and to be a maternally provided gene product. [provided by RefSeq, Apr 2009]