

Product datasheet for RC216642

TNNT3 (NM 006757) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: TNNT3 (NM_006757) Human Tagged ORF Clone

Tag: Myc-DDK TNNT3

Synonyms: beta-TnTF; DA2B2; TNTF

Mammalian Cell Neomycin

Selection:

E. coli Selection:

Symbol:

Vector: pCMV6-Entry (PS100001)

ORF Nucleotide >RC216642 representing NM_006757

Red=Cloning site Blue=ORF Green=Tags(s) Sequence:

Kanamycin (25 ug/mL)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC

GCCGCGATCGCC

ATGTCTGACGAGGAAGTTGAACAGGTGGAGGAGCAGTACGAAGAAGAAGAAGAAGCCCAGGAGGAAGAAGAGAG AAGTTCAAGAAGACACCGCAGAGGAGGACGCGGAAGAGGAGAAACCGAGACCCAAACTCACTGCTCCTAA GATCCCAGAAGGGGAGAAAGTGGACTTCGATGACATCCAGAAGAAGCGTCAGAACAAAGACCTAATGGAG CTCCAGGCCCTCATCGACAGCCACTTTGAAGCCCGGAAGAAGGAGGAGGAGGAGGAGCTGGTCGCTCTCAAAG AGAGAATCGAGAAGCGCCGTGCAGAGAGAGCGGAGCAGCAGAGGATTCGTGCAGAGAAGGAGAGGGAGCG CCAGAACAGACTGGCGGAGGAAAAGGCCAGAAGGGAGGAGGAGGATGCCAAGAGGAGGACGACGAC CTGAAGAAGAAGAAGCTCTGTCTTCCATGGGAGCCAACTACAGCAGCTACCTGGCCAAGGCTGACCAGA AGAGAGCAAGAAGCAGACAGCCCGGGAAATGAAGAAGAAGATTCTGGCTGAGAGACGCAAGCCGCTCAA CATCGATCACCTTGGTGAAGACAAACTGAGGGACAAGGCCAAGGAGCTCTGGGAGACCCTGCACCAGCTG GAGATTGACAAGTTCGAGTTTGGGGAGAAGCTGAAACGCCAGAAATATGACATCACCACGCTCAGGAGCC **GAAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT ACAAGGATGACGACGATAAGGTTTAA



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

TNNT3 (NM_006757) Human Tagged ORF Clone - RC216642

Protein Sequence: >RC216642 representing NM_006757

Red=Cloning site Green=Tags(s)

MSDEEVEQVEEQYEEEEEAQEEEEVQEDTAEEDAEEEKPRPKLTAPKIPEGEKVDFDDIQKKRQNKDLME LQALIDSHFEARKKEEEELVALKERIEKRRAERAEQQRIRAEKERERQNRLAEEKARREEEDAKRRAEDD LKKKKALSSMGANYSSYLAKADQKRGKKQTAREMKKKILAERRKPLNIDHLGEDKLRDKAKELWETLHQL EIDKFEFGEKLKRQKYDITTLRSRIDQAQKHSKKAGTPAKGKVGGRWK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja2280 f11.zip

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_006757

ORF Size: 774 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).





Reconstitution Method:

- 1. Centrifuge at 5,000xg for 5min.
- 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
- 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: <u>NM 006757.4</u>

RefSeq Size: 1217 bp RefSeq ORF: 777 bp


Locus ID: 7140
UniProt ID: P45378
Cytogenetics: 11p15.5
Domains: Troponin
MW: 31 kDa

Gene Summary: The binding of Ca(2+) to the trimeric troponin complex initiates the process of muscle


contraction. Increased Ca(2+) concentrations produce a conformational change in the troponin complex that is transmitted to tropomyosin dimers situated along actin filaments. The altered conformation permits increased interaction between a myosin head and an actin filament which, ultimately, produces a muscle contraction. The troponin complex has protein subunits C, I, and T. Subunit C binds Ca(2+) and subunit I binds to actin and inhibits actin-myosin interaction. Subunit T binds the troponin complex to the tropomyosin complex and is also required for Ca(2+)-mediated activation of actomyosin ATPase activity. There are 3 different troponin T genes that encode tissue-specific isoforms of subunit T for fast skeletal-, slow skeletal-, and cardiac-muscle. This gene encodes fast skeletal troponin T protein; also known as troponin T type 3. Alternative splicing results in multiple transcript variants encoding additional distinct troponin T type 3 isoforms. A developmentally regulated switch between fetal/neonatal and adult troponin T type 3 isoforms occurs. Additional splice variants have been described but their biological validity has not been established. Mutations in this gene may cause distal arthrogryposis multiplex congenita type 2B (DA2B). [provided by



Product images:




Circular map for RC216642



Western blot validation of overexpression lysate (Cat# [LY402019]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC216642 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).





Coomassie blue staining of purified TNNT3 protein (Cat# [TP316642]). The protein was produced from HEK293T cells transfected with TNNT3 cDNA clone (Cat# RC216642) using MegaTran 2.0 (Cat# [TT210002]).