

Product datasheet for **RC216662**

OOEP (NM_001080507) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: OOEP (NM_001080507) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: OOEP
Synonyms: C6orf156; FLOPED; HOEP19; KHDC2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC216662 representing NM_001080507
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTCGATGATGCTGGTGCCGCTGAGTCCCAGCGGGCAAACAGACTCCGGCCACTCCCTGGAGCAGC
TGCGTAGGTTACCACTTCCGCCGCCACAGATTCGCATCCGGCCCTGGTGGTTTCCGGTGCAGGAAGT
AGACCCCTTTGGTGTCTACCTAGAGGCATGGCTGGCAGACGAGCTCTTTGGCCAGACCGAGCCATAATT
CCAGAAATGGAGTGGACGAGCCAGGCCCTGCTGACAGTGGACATAGTGGACTCAGGGAACCTAGTCGAAA
TCACCGTTTTTCGGGCGGCCCGGTACAGAATCGGGTGAAGAGCATGCTCCTGTGCCTGGCATGGTTTCA
CCGAGAACATCGTGCCCGAGCTGAGAAGATGAAACACCTTGAGAAGAAGTGAAGGCCCATGCATCAGAC
CCCCACTCTCCCAGGATCCTGTTGCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC216662 representing NM_001080507
Red=Cloning site Green=Tags(s)

MVDDAGAAESQRGKQTPAHSLEQLRRLPLPPPQIRIRPWWFPVQELRDPLVFYLEAWLADELFGPDRAII
PEMEWTSQALLTVDIVDSGNLVEITVFGPRPVQNRVKSMLLCLAWFHREHRARAEMKHKLEKNLKAHASD
PHSPQDPVA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mg3696_d07.zip



[View online >](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001080507

ORF Size: 447 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001080507.2](#), [NP_001073976.1](#)

RefSeq Size: 681 bp


RefSeq ORF: 450 bp

Locus ID: 441161


UniProt ID: [A6NGQ2](#)
 Cytogenetics: 6q13
 MW: 17 kDa

Gene Summary: As a member of the subcortical maternal complex (SCMC), plays an essential role for zygotes to progress beyond the first embryonic cell divisions.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC216662



Western blot validation of overexpression lysate (Cat# [LY421062]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC216662 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).