

Product datasheet for **RC216668**

ALDH1A2 (NM_170697) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ALDH1A2 (NM_170697) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ALDH1A2
Synonyms:	RALDH(II); RALDH2; RALDH2-T
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC216668 representing NM_170697
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGATGCTTCAGAAAGGGGACGTCTGTTGGATAAGCTTGCAGACTTGGTGAACGGGACAGGGCAGTTC
 TTGCAACCATGGAATCCCTAAATGGTGGCAAACCATTCCTGCAAGCTTTTATGTGGATTTGCAGGGCGT
 CATCAAAACCTTTTCGATATTACGCAGGCTGGGCTGATAAAATTCATGGGATGACCATTCTGTAGATGGA
 GACTATTTTACCTTTACAAGACATGAACCCATTGGAGTGTGTGGACAGATCATCCATGGAACCTCCCCC
 TGCTGATGTTTGCCTGGAATAAGCTCCAGCTTTGTGCTGTGGCAATACAGTAGTTATTAAGCCAGCAGA
 GCAAACACCCTCAGTGCCTCTACATGGGAGCCCTCATCAAGGAGGCTGGCTTTCCTCCGGGGTCATC
 AATATTTTCCAGGATATGGCCAACGGCTGGGCGAGCAATAGCTTCTCACATTGGCATAGACAAGATTG
 CATTACAGGGTCTACTGAGTTGGAAAGCTTATCCAAGAAGCAGCTGGAAGAAGTAATTTGAAGAGAGT
 AACTCTGGAACCTGGAGGCAAAGTCTAATATATTTTTGCTGATGCTGACTTGGACTATGCTGTGGAG
 CAGGCCACCAGGGTGTGTTCTTCAATCAAGGTCAGTGCCTGCAGGCTCTCGCATCTTCGTGGAGG
 AGTCCATCTATGAGGAGTTTGTGAGAAGAAGCGTGGAGCGGGCCAAGAGGCGCGTGTGGGGAGTCCCTT
 TGACCCACCCTGAGCAGGGTCCCGAGATTGATAAGAAACAGTACAACAAGATCTTGGAACTCATCCAG
 AGTGGTGTGGCTGAGGGCGCCAAGCTGGAATGTGGAGGCAAAGGACTGGGCCGAAAGGGTTTTTCATTG
 AGCCACAGTGTTCACACGTCCTGATGATATGCGGATTGCCAAGGAGGAGATCTTTGGCCCTGTTCA
 GGAAATTTGAGATTTAAGACGATGGATGAAGTTATCGAAAGAGCCAATAACTCAGACTTTGGACTCGTA
 GCAGCTGTCTTACTAATGACATCAACAAGGCCCTCACAGTGTCTTCTGCAATGCAAGCTGGGACTGTTT
 GGATCAATTGTTACAATGCCTTAAATGCCAGAGCCCTTTGGGGATTCAAGATGTCTGAAATGGGAG
 AGAAATGGGAGAATTTGGCTTCCGGGAGTACTCAGAAGTTAAGACGGTGACAGTAAAGATCCCCAGAAG
 AACTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC216668 representing NM_170697
 Red=Cloning site Green=Tags(s)

MDASERGRLLDKLADLVERDRAVLATMESLNGGKPFLLQAFYVDLQGVIKTFRYYAGWADKIHGMTIPVDG
 DYFTFTRHEPIGVCGQIIPWNFPLLMFAWKIAPALCCGNTVVIKPAEQTPLSALYMGALIKEAGFPPGVI
 NILPGYGPTAGAAIASHIGIDKIAFTGSTEVGKLIQEAAGRSNLKRVTLLEGGKSPNIIFADADLDYAVE
 QAHQGVFFNQGCCTAGSRIFVEESIYEEFVRRSVERAKRRVVGSPFPDPTTEQGPQIDKKQYNKILELIQ
 SGVAEGAKLECGKGLGRKGFIEPTVFSNVTDDMRIAKEEIFGPVQEILRFKTMDEVIERANNSDFGLV
 AAVFTNDINKALTVSSAMQAGTVWINCYNALNAQSPFGGFKMSGNGREMGEGFLREYSEVKTVTVKIPQK
 NS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/mk6595_c10.zip

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:



ACCN: NM_170697

ORF Size: 1266 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_170697.3](#)

RefSeq Size: 3142 bp

RefSeq ORF: 1269 bp

Locus ID: 8854

UniProt ID: [O94788](#)

Cytogenetics: 15q21.3

Protein Families: Druggable Genome


Protein Pathways: Metabolic pathways, Retinol metabolism

MW: 46 kDa


Gene Summary:

This protein belongs to the aldehyde dehydrogenase family of proteins. The product of this gene is an enzyme that catalyzes the synthesis of retinoic acid (RA) from retinaldehyde. Retinoic acid, the active derivative of vitamin A (retinol), is a hormonal signaling molecule that functions in developing and adult tissues. The studies of a similar mouse gene suggest that this enzyme and the cytochrome CYP26A1, concurrently establish local embryonic retinoic acid levels which facilitate posterior organ development and prevent spina bifida. Four transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, May 2011]

Product images:



Circular map for RC216668



Western blot validation of overexpression lysate (Cat# [LY406883]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC216668 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).