

Product datasheet for **RC216669L4V**

DISC1 (NM_001012957) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DISC1 (NM_001012957) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DISC1
Synonyms:	C1orf136; SCZD9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001012957
ORF Size:	2496 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216669).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001012957.1 , NP_001012975.1
RefSeq Size:	7003 bp
RefSeq ORF:	2499 bp
Locus ID:	27185
UniProt ID:	Q9NRI5
Cytogenetics:	1q42.2
MW:	90.9 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein with multiple coiled coil motifs which is located in the nucleus, cytoplasm and mitochondria. The protein is involved in neurite outgrowth and cortical development through its interaction with other proteins. This gene is disrupted in a t(1;11) (q42.1;q14.3) translocation which segregates with schizophrenia and related psychiatric disorders in a large Scottish family. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]