

Product datasheet for RC216679

OR2F2 (NM_001004685) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OR2F2 (NM_001004685) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	OR2F2
Synonyms:	OR7-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC216679 representing NM_001004685 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAATAGATAAACCAGACGTGGGTGAGAGAATTTATTCTCCTTGGCTTATCCAGTGACTGGTGCCTC
AGATATCCCTGTTTTCCCTGTTCTTGGTCACATACCTCATGACAGTGCTGGGAACTGTCTCATTGTCTC
TCTGATCAGACTGGACAGCCGACTCCACACTCCCATGATTTCTTCTCACCAACCTCTCCCTTGTGAT
GTCTCCTATGCCACAAGCGTAGTCCCCAGCTGCTGGCACATTTCTTGCAGAACATAAAGCCATCCCAT
TCCAGAGCTGTGCAGCCAGTTATTTTTCTCCCTGGCCTTGGGTGGGATTGAGTTTGTCTCCTGGCAGT
GATGGCCTATGACCGCCATGTGGCTGTGTCTGACCGCCTGCGATACTCGGCCATCATGCATGGAGGGCTG
TGTGCTAGGTTGGCCATCACATCCTGGTCACTGGCTCCATCAACTCTCTTGTGCAGACTGCTATACCT
TTCAGTGTCCCATGTGCACTAACAAGTTTATTGATCACATATCCTGTGAACCTCTAGCTGTGGTCAGGCT
GGCTTGTGTGGACACCTCCTCCAATGAGGCTGCCATCATGGTGTCTAGCATTGTTCTTCTGATGACACCT
TTCTGCCTGGTTCTGTTGTCTACATCCGGATCATCTCCACCATCCTAAAGATCCAGTCCAGAGAAGGAA
GAAAGAAAGCCTCCACACGTGTGCCTCTACCTCACGGTGGTGGCCTGTGCTACGGCACAAACGATTTT
CACTTACATCCAGCCCCACTCTGGTCCCTCAGTCTTCAAGAGAAGCTGATCTCTGTCTTCTATGCCATT
GTTATGCCTCTGCTGAACCCTGTGATTTATAGTCTAAGGAATAAAGAGGTGAAGGGGCCCTGGCATAAAC
TATTAGAGAAATCTCTGGGTTAACATCCAAGCTGGGAACT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC216679 representing NM_001004685
Red=Cloning site Green=Tags(s)

```
MEIDNQTWVREFILLGLSSDWCTQISLFSFLVTVLMTVLGNCLIVLLIRLDSRLHTPMYFFLTNLSLVD
VSYATSVVPQLLAHFLAEHKAIPFQSCAAQLFFSLALGGIEFVLLAVMAYDRHVAVSDRLRYSAIMHGGL
CARLAITSWVSGSINSLVQTAITFQLPMCTNKFIDHISCELLAVVRLACVDTSSNEAAIMVSSIVLLMTP
FCLVLLSYIRIIISTILKIQSREGRKKAFTCASHLTVVVALCYGTTIFTYIQPHSGPSVLQEKLISVFYAI
VMPLLNPIYSLRNKEVKGAWHKLLEKFSGLTSKLG
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8002_d02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001004685

ORF Size: 951 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001004685.1 , NP_001004685.1
RefSeq Size:	954 bp
RefSeq ORF:	954 bp
Locus ID:	135948
UniProt ID:	O95006
Cytogenetics:	7q35
Protein Families:	Transmembrane
Protein Pathways:	Olfactory transduction
MW:	35.1 kDa
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC216679