

Product datasheet for **RC216681L3V**

ARHGEF9 (NM_015185) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ARHGEF9 (NM_015185) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ARHGEF9
Synonyms:	COLLYBISTIN; DEE8; EIEE8; HPEM-2; PEM-2; PEM2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_015185
ORF Size:	1548 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216681).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015185.1
RefSeq Size:	5413 bp
RefSeq ORF:	1551 bp
Locus ID:	23229
UniProt ID:	O43307
Cytogenetics:	Xq11.1
Domains:	RhoGEF, SH3, PH
MW:	60.8 kDa


[View online »](#)

Gene Summary:

The protein encoded by this gene is a Rho-like GTPase that switches between the active (GTP-bound) state and inactive (GDP-bound) state to regulate CDC42 and other genes. This brain-specific protein also acts as an adaptor protein for the recruitment of gephyrin and together these proteins facilitate receptor recruitment in GABAergic and glycinergic synapses. Defects in this gene are the cause of startle disease with epilepsy (STHEE), also known as hyperekplexia with epilepsy, as well as several other types of cognitive disability. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2017]