

Product datasheet for **RC216706A1V**

Human ICAM2 (NM_001099788) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ICAM2 (NM_001099788) AAV Particle
Tag:	Myc-DDK
Symbol:	ICAM2
Synonyms:	CD102
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC216706 representing NM_001099788 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCCTCTTTCGGTTACAGGACCCTGACTGTGGCCCTCTTACCCTGATCTGCTGTCCAGGATCGGATG
AGAAGGTATTCGAGGTACACGTGAGGCCAAAGAAGCTGGCGGTTGAGCCAAAGGGTCCCTCGAGGTCAA
CTGCAGCACCACCTGTAACAGCCTGAAGTGGTGGTCTGGAGACCTCTCTAGATAAGATTCTGCTGGAC
GAACAGGCTCAGTGAAACATTACTTGGTCTCAAACATCTCCCATGACACGGTCCCTCAATGCCACTTCA
CCTGCTCCGGAAGCAGGAGTCAATGAATTCCAACGTCAGCGTGTACCAGCCTCCAAGGCAGGTATCCTCT
GACACTGCAACCACTTTGGTGGCTGTGGCAAGTCCTTACCATTGAGTGCAGGGTGCCACCGTGGAG
CCCCTGGACAGCCTCACCTCTTCTGTTCGTTGGCAATGAGACTCTGCACTATGAGACCTTCGGAAGG
CAGCCCCTGCTCCGAGGAGGCCACAGCCACATTCAACAGCACGGCTGACAGAGAGGATGGCCACCGCAA
CTTCTCTGCCTGGCTGTGCTGGACTTGATGTCTCGCGGTGGCAACATCTTTCACAAACACTCAGCCCCG
AAGATGTTGGAGATCTATGAGCCTGTGTCGGACAGCCAGATGGTCATCATAGTCACGGTGGTGTGCGGTG
TGCTGTCCCTGTTCTGTGACATCTGTCTGCTCTGCTTCATCTTCGGCCAGCACTTGCGCCAGCAGCGGAT
GGGCACCTACGGGGTGCAGCGGCTTGAGGAGGCTGCCCCAGGCCTTCGGCCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC216706 representing NM_001099788 Red=Cloning site Green=Tags(s) MSSFGYRTLTVLFTLICCPGSDEKVFVHVRPKKLAVEPKGSLEVNCSTTCNQPEVGGLTSLDKILLD EQAQWKHYLVSNISHDTVLQCHFTCSGKQESMNSNVSVYQPPRQVILTLQPTLVAVGKSFTIECRVPTVE PLDSLTLFLFRGNETLHYETFGKAAPAPQEATATFNSTADREDGHRNFSCLAVLDLMSRGGNIFHKHSAP KMLEIYEPVSDSQMVIIIVTVSVLLSLFVTSVLLCFIFGQHLRQQRMGTYGVRRAAWRRLPQAFRP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001099788
ORF Size:	825 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001099788.1, NP_001093258.1</u>
RefSeq Size:	1235 bp
RefSeq ORF:	828 bp
Locus ID:	3384
UniProt ID:	<u>P13598</u>
Cytogenetics:	17q23.3
MW:	30.7 kDa