

## Product datasheet for **RC216731L3V**

### LARP4 (NM\_052879) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	LARP4 (NM_052879) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LARP4
Synonyms:	PP13296
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_052879
ORF Size:	2172 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216731).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_052879.4</a>
RefSeq Size:	6537 bp
RefSeq ORF:	2175 bp
Locus ID:	113251
UniProt ID:	<a href="#">Q71RC2</a>
Cytogenetics:	12q13.12
Protein Families:	Transcription Factors
MW:	80.6 kDa



[View online »](#)

**Gene Summary:**

RNA binding protein that binds to the poly-A tract of mRNA molecules (PubMed:21098120). Associates with the 40S ribosomal subunit and with polysomes (PubMed:21098120). Plays a role in the regulation of mRNA translation (PubMed:21098120). Plays a role in the regulation of cell morphology and cytoskeletal organization (PubMed:21834987, PubMed:27615744). [UniProtKB/Swiss-Prot Function]