

Product datasheet for **RC216738A1V**

Human CD42c (GP1BB) (NM_000407) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CD42c (GP1BB) (NM_000407) AAV Particle
Tag: Myc-DDK
Symbol: CD42c
Synonyms: BDPLT1; BS; CD42C; GPIBB; GPIbbeta
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC216738 representing NM_000407
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGCTCCGGGCGCGCGGGGCGCTGAGCTTACTGCTCCTGCTGCTGGCCCCGCCGAGCCGCCCGGCCG
CAGGTTGCCGGCGCCCTGTAGCTGCGCGGGGACGCTCGTGGACTGCGGGCGCCGCGGGCTGACTTGGGC
CTCGCTGCCGACCGCTTCCCTGTGACACAACCGAGCTGGTGTGACCGGCAACAACCTGACGGCGCTG
CCGCCGGGGCTGCTGGACGCGCTGCCGCGCTGCGCACCGCACACCTGGGCGCCAACCCCTGGCGCTGCG
ACTGCCGCTTGTGCGCTGCGCGCTGGTGGCGGGCGCCCCGAGCGTGCCTACCGCGACCTGCG
TTGCGTGGCGCCCCAGCGCTGCGCGGGCGCTGCTGCCCTATCTGGCCGAGGACGAGCTGCGCGCGCT
TGCGCTCCCGGGCGCTGCTGGGGGCGCTGGCGGCGCAGCTTGCCTGCTGGGCTTGGGCTGCTGC
ACGCGTTGCTGCTGGTGTGCTGTGCGCGCTGCGGAGGCTGCGGGCCCGGGCCCGCTGCGCGCGC
AGCCCGGCTGTCGCTGACCGACCGCTGGTGGCCGAGCGAGCCGAACCGACGAGTCC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC216738 representing NM_000407
Red=Cloning site Green=Tags(s)

MGSGPRGALSLLLLLAPPSPRAAGCPAPCSCAGTLVDCGRRGLTWASLPTAFPVDTTELVL TGNL TAL
PPGLLDALPALRTAHLGANPWRCDCRLVPLRAWLAGRPERAPYRDLRCVAPPALRGRLLPYLADELRAA
CAPGPLCWGALAAQLALLGLLHALLLVLLCRLRRLRARARARAAARLSLTDPLVAERAGTDES

SGPTRTRRLE**QKLI**SEEDLAANDILDYKDDDDKV

Species: Human



Serotype:	AAV-2
ACCN:	NM_000407
ORF Size:	618 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_000407.3
RefSeq Size:	955 bp
RefSeq ORF:	621 bp
Locus ID:	2812
UniProt ID:	P13224
Cytogenetics:	22q11.21
MW:	21.72 kDa