

Product datasheet for **RC216741L3V**

GAMT (NM_138924) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GAMT (NM_138924) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GAMT
Synonyms:	CCDS2; HEL-S-20; PIG2; TP53I2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_138924
ORF Size:	807 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216741).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138924.2
RefSeq Size:	1787 bp
RefSeq ORF:	810 bp
Locus ID:	2593
UniProt ID:	Q14353
Cytogenetics:	19p13.3
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways:	Arginine and proline metabolism, Glycine, serine and threonine metabolism, Metabolic pathways
MW:	29.8 kDa
Gene Summary:	<p>The protein encoded by this gene is a methyltransferase that converts guanidoacetate to creatine, using S-adenosylmethionine as the methyl donor. Defects in this gene have been implicated in neurologic syndromes and muscular hypotonia, probably due to creatine deficiency and accumulation of guanidinoacetate in the brain of affected individuals. Two transcript variants encoding different isoforms have been described for this gene. Pseudogenes of this gene are found on chromosomes 2 and 13. [provided by RefSeq, Feb 2012]</p>