

## Product datasheet for **RC216749**

### **PDE9A (NM\_001001567) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	PDE9A (NM_001001567) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PDE9A
Synonyms:	HSPDE9A2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>RC216749 representing NM\_001001567  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGGATCCGGCTCCTCCAGCTACCGGCCAAGGCCATCTACCTGGACATCGATGGACGCATTGAGAAGG  
 TAATCTTCAGCAAGTACTGCAACTCCAGCGACATCATGGACCTGTTCTGCATCGCCACCGCCTGCCTCG  
 GAACACGACCATCTCCCTGCTGACCACCGACGACGCCATGGTCTCCATCGACCCACCATGCCCGGAAT  
 TCAGAACGCACTCCGTACAAAGTGAACCTGTGGCCATCAAGCAACTCTCCGAGAGAGAAGAATTAATCC  
 AGAGCGTGTGGCGCAGTTGACAGCAGTTCTCAAGAGCATTAAAAATCAATGAACTGAAAGCTGAAGT  
 TGCAAATCACTTGGCTGTCTAGAGAAACGCGTGAATTTGGAAGGACTAAAAGTGGTGGAGATTGAGAAA  
 TGCAAGAGTGACATTAAGAAGATGAGGGAGGAGCTGGCGCCAGAAGCAGCAGGACCAACTGCCCTGTA  
 AGTACAGTTTTTTGGATAACCACAAGAAGTTGACTCCTCGACCGATGTTCCCACTTACCCCAAGTACCT  
 GCTCTCTCCAGAGACCATCGAGGCCCTGCGGAAGCCGACCTTTGACGTCTGGCTTTGGGAGCCCAATGAG  
 ATGCTGAGCTGCCTGGAGCACATGTACCACGACCTCGGGCTGGTCAAGGACTTCAGCATCAACCCTGTCA  
 CCCTCAGGAGGTGGCTGTTCTGCGTCCACGACAACACGAAACAAACCCCTTCCACAACCTCCGGCACTG  
 CTTCTGCGTGGCCAGATGATGTACAGCATGGTCTGGCTCTGCAGTCTCCAGGAGAAGTTCTCACAACG  
 GATATCCTGATCCTAATGACAGCGGCCATCTGCCACGATCTGGACCATCCCGGCTACAACAACACGTACC  
 AGATCAATGCCCGCACAGAGCTGGCGGTCCGCTACAATGACATCTCACCGCTGGAGAACCACCACTGCGC  
 CGTGGCCTTCCAGATCCTCGCCGAGCCTGAGTGAACATCTTCTCCAACATCCCACCTGATGGGTTCAAG  
 CAGATCCGACAGGAATGATCACATTAATCTTGGCCACTGACATGGCAAGACATGCAGAAATTAATGGATT  
 CTTTCAAAGAGAAAAATGGAGAATTTTACTACAGCAACGAGGAGCACATGACCCTGCTGAAGATGTTTT  
 GATAAAATGCTGTGATATCTTAACGAGGTCCGTCCAATGGAAGTCGAGAGCCTTGGGTGGACTGTTTA  
 TTAGAGGAATATTTTATGCAGAGCGACCGTGAGAAGTCAGAAGGCCTTCTGTGGCACCGTTTATGGACC  
 GAGACAAAGTGACCAAGGCCACAGCCAGATTGGGTTTCAAGTTTGTCTGATCCCAATGTTTGAAC  
 AGTGACCAAGCTTCCCATGGTTGAGGAGATCATGCTGCAGCCACTTTGGGAATCCCAGATCGCTAC  
 GAGGAGCTGAAGCGGATAGATGACGCCATGAAAGATTACAGAAGAAGACTGACAGCTTGACGTCTGGGG  
 CCACCGAGAAGTCCAGAGAGAGAAGCAGAGATGTAAAAACAGTGAAGGAGACTGTGCC

**ACGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGA**  
**TTACAAGGATGACGACGATAAGGTTAA**

**Protein Sequence:**

>RC216749 representing NM\_001001567  
 Red=Cloning site Green=Tags(s)


MGGSSSYRPKAIYLDIDGRIQKVIKSKYCNSSDIMDLFCIATGLPRNTTISLLTTDDAMVSIPTMPAN  
 SERTPYKVRPVAIKQLSEREELIQSVLAQVAEQFSRAFKINELKAEVANHLAVLEKRVELEGLKVVEIEK  
 CKSDIKMREELAARSSRTNCPCKYSFLDNHKKLTPRRDVPTYPKYLLSPETIEALRKPTFDVWLWEPNE  
 MLSCLEHMYHDLGLVRDFSINPVTLLRRWLFVVDNYRNNPFHFRHCFVQMMYSMVWLCSEKFSQT  
 DILILMTAAICHDLDPGYNNTYQINARTELAVRYNDISPLENHHCVAFAQILAEPECNIFSNIPPDGFK  
 QIRQGMITLILATDMARHAEIMDSFKEKMFNFYSNEEHMTLLKMILIKCCDISNEVRPMEVAEPWVDC  
 LEEYFMQSDREKSEGLPVAPFMDRDKVTKATAQIGFIKFLIPMFETVTKLFPMVVEIMLQPLWESRDY  
 EELKRIDDAMKELQKKTDSLTSGATEKSRERSRDVKNSEGDC

**TRRLEQKLISEEDLAANDILDYKDDDDKV**


**Restriction Sites:**

Sgfl-NotI

Cloning Scheme:



Plasmid Map:



<b>ACCN:</b>	NM_001001567
<b>ORF Size:</b>	1599 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001001567.1</a> , <a href="#">NP_001001567.1</a>
<b>RefSeq Size:</b>	1923 bp
<b>RefSeq ORF:</b>	1602 bp
<b>Locus ID:</b>	5152
<b>UniProt ID:</b>	<a href="#">O76083</a>
<b>Cytogenetics:</b>	21q22.3
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Progesterone-mediated oocyte maturation, Purine metabolism
<b>MW:</b>	61.5 kDa
<b>Gene Summary:</b>	The protein encoded by this gene catalyzes the hydrolysis of cAMP and cGMP to their corresponding monophosphates. The encoded protein plays a role in signal transduction by regulating the intracellular concentration of these cyclic nucleotides. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]