

Product datasheet for **RC216763L3V**

Glycoprotein 2 (GP2) (NM_001502) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Glycoprotein 2 (GP2) (NM_001502) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | Glycoprotein 2 |
| Synonyms: | ZAP75 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001502 |
| ORF Size: | 1602 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC216763). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001502.2 |
| RefSeq Size: | 2439 bp |
| RefSeq ORF: | 1605 bp |
| Locus ID: | 2813 |
| UniProt ID: | P55259 |
| Cytogenetics: | 16p12.3 |
| Domains: | zona_pellucida |
| Protein Families: | Druggable Genome, Secreted Protein, Transmembrane |



[View online »](#)

MW: 58.9 kDa

Gene Summary: This gene encodes an integral membrane protein that is secreted from intracellular zymogen granules and associates with the plasma membrane via glycosylphosphatidylinositol (GPI) linkage. The encoded protein binds pathogens such as enterobacteria, thereby playing an important role in the innate immune response. The C-terminus of this protein is related to the C-terminus of the protein encoded by the neighboring gene, uromodulin (UMOD). Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2015]