

## Product datasheet for **RC216777L4V**

### EPB41L3 (NM\_012307) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	EPB41L3 (NM_012307) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EPB41L3
Synonyms:	4.1B; DAL-1; DAL1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012307
ORF Size:	3261 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216777).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012307.2</a> , <a href="#">NP_036439.2</a>
RefSeq Size:	4446 bp
RefSeq ORF:	3264 bp
Locus ID:	23136
UniProt ID:	<a href="#">Q9Y2J2</a>
Cytogenetics:	18p11.31
Domains:	B41, SAB
Protein Pathways:	Tight junction



[View online »](#)

**MW:** 120.5 kDa

**Gene Summary:** Tumor suppressor that inhibits cell proliferation and promotes apoptosis. Modulates the activity of protein arginine N-methyltransferases, including PRMT3 and PRMT5.  
[UniProtKB/Swiss-Prot Function]