

## Product datasheet for **RC216813L2V**

### Oct-2 (POU2F2) (NM\_002698) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Oct-2 (POU2F2) (NM_002698) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Oct-2
Synonyms:	Oct-2; OCT2; OTF2
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_002698
ORF Size:	1389 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC216813).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_002698.1</a>
RefSeq Size:	2048 bp
RefSeq ORF:	1392 bp
Locus ID:	5452
UniProt ID:	<a href="#">P09086</a>
Cytogenetics:	19q13.2
Domains:	POU, homeobox
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 49.3 kDa

**Gene Summary:** The protein encoded by this gene is a homeobox-containing transcription factor of the POU domain family. The encoded protein binds the octamer sequence 5'-ATTGCAT-3', a common transcription factor binding site in immunoglobulin gene promoters. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]